


Callisto 930 nm OCT Imaging System



Overview | Software | OCT Tutorial | Brochures | Feedback | Selection Guide

Features

- Configurable OCT Systems Optimized for High-Sensitivity Imaging
 - 1.2 kHz A-Scan Rate and 107 dB Sensitivity
 - 7 μ m Axial Resolution in Air with 1.7 mm Imaging Depth (930 nm Center Wavelength)
- Includes Computer and ThorImage[®] OCT Software Package (See the Software Tab)
- Build-Your-Own and Preconfigured Systems Available
- See the Brochures Tab for More Information on Thorlabs' OCT Systems

Choose Components to Build or Customize Your OCT System

- One 930 nm, High-Sensitivity Base Unit (Required)
- Standard, User-Customizable, and Handheld Scanners Available
- Scan Lens Kits to Optimize Lateral Resolution and Focal Length for Your Application
- Ring- and Immersion-Style Sample Z-Spacers for Air or Liquid Imaging Applications
- Scanner Stand and Translation Stage Accessories
- Contact Our [OCT Team](#) to Request a Quote and Discuss Building a System

Optical Coherence Tomography (OCT) is a noninvasive optical imaging technique that produces real-time, 2D cross-sectional and 3D volumetric images of a sample. This technique provides structural information about the sample based on light backscattered from different layers of material within that sample, producing images with micron-level resolution and millimeters of imaging depth. OCT imaging can be considered as an optical analog to ultrasound imaging that achieves higher resolution at the cost of decreased penetration depth. In addition to high resolution, the non-contact, noninvasive nature of OCT makes it well suited for imaging samples such as biological tissue, small animals, and industrial materials.

Thorlabs' Callisto Series of OCT Imaging Systems are designed for high sensitivity and are ideal for imaging static, *in vitro* samples. The 64-bit software pre-installed on the included computer displays and processes 2D and 3D OCT data in real time. Choose from a number of scanner options including a robust standard scanner, user-customizable scanner, and the portable handheld scanner. Optional accessories are available below to customize your OCT system to meet the requirements of your application. Additionally, Thorlabs offers a complete, preconfigured OCT system for 930 nm based on the components sold on this page.

The components below can also be used to upgrade your existing Thorlabs OCT system with additional features and are fully compatible out of the box with Thorlabs' OCT systems and accessories. While most systems are upgradeable, we recommend contacting the [OCT Team](#) to determine the optimal solution for your system and intended application.



OCT Applications Team Based in Lübeck, Germany

We're Happy to Assist!

[Contact Us](#)

Exploring the Options?

We can provide recommendations based on your needs and partner with you to obtain images of samples provided by you demonstrating the effects of various components on image quality. Demos of our OCT systems can be arranged at our Sterling, VA (USA); Shanghai, China; Tokyo, Japan; and Lübeck, Germany facilities.

In the Budgetary Phase?

System prices vary based on the exact components. Through our conversations, we can ensure your system quote is tailored to your requirements.

OEM or Custom Projects?

Click [here](#) to learn about our OEM capabilities.

OCT System Comparison



Customization Guide



Application Articles

Click on the Image Below or in the Table to the Right for Details on Customization Options



Components and Accessories

- [OCT Base Unit](#) (Computer Included)
- [Scanning System](#)
- [Scan Lens Kit](#)
- [Reference Length Adapter](#) (For Standard Scanners Only)
- [Sample Z-Spacer](#)
- [Adjustable Scanner Stand](#)
- [Translation Stage](#)
- [Preconfigured System](#) (Z-Spacer Not Included)

Callisto Series Complete Preconfigured System



- Complete Preconfigured 930 nm OCT Imaging System
- High-Sensitivity Imaging (See Tables Below)
- Fully Customizable Using Other Callisto Series Components

Thorlabs offers this complete, preconfigured Callisto OCT system, which is fully compatible with all Callisto Series OCT components. The CAL110C1 configuration features a center wavelength of 930 nm and is designed for high-sensitivity imaging applications.

This Callisto Series OCT system configuration is built completely from components sold in sections located lower on this page. This preconfigured system includes the three mandatory OCT system core components (the base unit, a scanning system with its reference length adapter, and a scan lens kit), as well as two optional accessories (scanner stand and translation stage). For more information about a component included in the preconfigured system, click on the component description in the table to the lower left to navigate down to the related section on this page.

For information about the CAL110C1 system or to inquire about custom configurations, please contact oct@thorlabs.com.

Preconfigured System Included Components	
System Item #	CAL110C1
Base Unit	CAL110
Scanning System	OCTG-900 (Standard Scanner)
Scan Lens Kit	OCT-LK3-BB
Reference Length Adapter	OCT-RA3
Accessories: Stand and Stage	OCT-STAND(M) (Scanner Stand) and OCT-XYR1(M) (Translation Stage)

a Click on the component description to navigate down to the related section on this page.

Preconfigured System Key Specifications	
System Item #	CAL110C1
Imaging Depth (Air/Water)	1.7 mm / 1.3 mm
Axial Resolution (Air/Water)	7.0 μ m / 5.3 μ m
Lateral Resolution	8 μ m
A-Scan/Line Rate	1.2 kHz
Sensitivity	107 dB

Based on your currency / country selection, your order will ship from Newton, New Jersey

+1 Qty	Docs	Part Number - Universal	Price	Available
<input type="text" value="1"/>		CAL110C1 Spectral Domain OCT System, 930 nm, 7.0 μ m Resolution, 1.2 kHz	\$37,131.50	Lead Time
Add To Cart				

OCT Base Unit (Required OCT System Component)



- 930 nm High-Sensitivity Base Unit
- 107 dB Sensitivity at 1.2 kHz A-Scan Rate
- 1.7 mm Imaging Depth and 7 μ m Axial Resolution

The imaging performance of any OCT system is largely dependent on the design and components incorporated into the base unit. All of Thorlabs' OCT Base Units include an OCT engine, high-performance computer, pre-installed software, and a software development kit (SDK). For the Callisto OCT Base Unit, the engine is comprised of a superluminescent diode light source, scanning electronics, and a linear CCD array-based spectrometer for detection. The engine and detection components are integrated into a 420 mm x 320 mm x 149 mm (16.5" x 12.6" x 5.86") unit. For a fully operational system, one [scanning system](#) and a [scan lens kit](#) must be purchased along with a base unit (sold separately below).

Thorlabs' CAL110 Base Unit is optimized for high sensitivity (107 dB) at an A-Scan line rate of 1.2 kHz. Thus, the Callisto is ideal for imaging static, *in vitro* samples that require good contrast to distinguish features. The CAL110 offers up to 1.7 mm imaging depth with 7 μ m of axial resolution.

Base Unit Item #	CAL110
Center Wavelength	930 nm
Imaging Depth (Air/Water)	1.7 mm / 1.3 mm
Axial Resolution (Air/Water)	7.0 μ m / 5.3 μ m
A-Scan Line Rate	1.2 kHz
Sensitivity ^a	107 dB (at 1.2 kHz)
Maximum Pixels per A-Scan	512
Compatible Scanners	OCTP-900(M), OCTG-900, and OCTG-900

a Typical Values Measured Using a Scanner with a Common Reference/Sample Path and 50% Path Split

**Also Required:
Scanning System
Scan Lens Kit**

To be functional, an OCT system build must include a base unit, a scanning system, and a scan lens kit.

Computer Specifications	
Operating System	Windows 7: 64 Bit
Processor	Dual Core >2.5 GHz
Memory	4 GB
Hard Drive	>250 GB
Data Acquisition	USB

a. Computer Specifications Subject to Change
b. We are in the process of switching the pre-installed operating system; some computers may come with Windows 10 installed upon delivery. This will not affect the performance of the OCT software.

Based on your currency / country selection, your order will ship from Newton, New Jersey

*1	Qty	Docs	Part Number - Universal		Price	Available
	<input type="text"/>		CAL-110	Callisto OCT Base Unit, 930 nm, 7.0 µm Resolution, 1.2 kHz	\$20,157.10	Lead Time

[Add To Cart](#)

Scanning Systems (Required OCT System Component)

- ▶ Scan an OCT Light Source Beam Across a Sample to Acquire 2D or 3D Images
- ▶ Three Available Options
 - Standard Scanner for High Stability and Ease-of-Use
 - User-Customizable Scanners with Open Construction for Customization of Scan Path
 - Portable Handheld Scanner for Applications Requiring High Mobility

Thorlabs' OCT Scanning Systems are designed to scan the OCT light source beam across a sample for 2D cross-sectional and 3D volumetric imaging. OCT applications can vary widely, from live animal imaging to industrial materials analysis, with each requiring a different set of scanning parameters. We currently offer three types of beam scanning systems for use with our Callisto Base Unit: standard, user-customizable, and handheld.

Each scanner contains an OCT interferometer with a sample arm and a reference arm. The reference arm of the OCT interferometer is placed near the sample and housed within the scanning system itself to guarantee the phase stability of the sample arm relative to the reference arm. To account for different sample distances and reflectivities (e.g., while imaging through water), the reference arm path length, as well as the reference arm intensity, is adjustable. To minimize image distortion caused by dispersion, our OCT systems are designed to optically match the reference and sample arm lengths to the greatest extent possible. Dispersion effects from the sample (e.g., imaging through water or glass) can be compensated for using the included ThorImage OCT software. For customers interested in dual-path setups, any of these scanners can be configured without a beamsplitter; please contact sales@thorlabs.com for more information.

All scanners are equipped with an integrated camera that can obtain real-time, en face-video of the sample during OCT measurements when used with our ThorImage OCT software (see the Software tab for details). Illumination of the sample is provided by a ring of user-adjustable white light LEDs around the exit aperture of each scanner.

Standard Scanner

The OCTG-900 Standard Scanner is ideal for imaging applications that require a stable, easy-to-operate setup. The entire design of the standard scanner is contained within a rugged, light-tight housing that minimizes the risk of misalignment. The standard scanner is equipped with a reference path length distance indicator for ease-of-use during reference adjustments. A knob located at the top of the standard scanner allows for fine adjustments to the reference path length.

User-Customizable Scanner

The OCTP-900(M) User-Customizable Scanner is designed with an open construction to enable easy customization of the optical beam path using Thorlabs' standard optomechanical components. This scanner features SM1 (1.035" x 40) ports and 4+40 tapped holes at several locations that allow mounting of SM1-threaded or 30 mm cage-compatible components, respectively. The scan lens port is directly compatible with either RMS 3 x 7.5 or SM1-threaded components, and can be converted to other thread standards, such as RMS (0.800"-36) via our selection of [thread adapters](#). Additional scanning and non-scanning optical input/output ports are available for integration of a laser for fluorescence excitation or additional sample illumination.

Compact Handheld Scanner

The compact and lightweight OCTH-900 Handheld Scanner is specifically designed for applications requiring high mobility. Easy access buttons located directly on the scanner enable fingertip control of our ThorImage OCT Software. Users can program each button from a selection of imaging and acquisition software controls and the software uses visual and audio feedback for button presses. The OCTH-900 features a removable cover under the handle that provides access to the reference length and intensity settings. Compatible scan lens kits and sample z-spacers for the OCTH-900 are sold below; z-spacers help maintain the correct working distance when using the handheld scanner. Please note that due to the limitations of the internal MEMS scanner, the frame rate (i.e., B-Scan rate) is limited to 25 frames per second when using the handheld scanner.



Scanner Type	Item #	Compatible Base Unit
Standard Scanner ^a	OCTG-900	CAL110
User-Customizable Scanner	OCTP-900(M)	
Handheld Scanner ^b	OCTH-900	

a. Standard Scanner Requires Purchase of [Reference Length Adapter](#)
b. [Sample Z-Spacers](#) Recommended When Using OCTH-900

**Also Required:
Base Unit
Scan Lens Kit**

To be functional, an OCT system build must include a base unit, a scanning system, and a scan lens kit.

[Click to Enlarge](#)
OCTH-900 Handheld Scanner with OCTH-900 Sample Z-Spacer

Based on your currency / country selection, your order will ship from Newton, New Jersey

*1	Qty	Docs	Part Number - Imperial		Price	Available
	<input type="text"/>		OCTP-900	Customer Inspired! User-Customizable Scanner for 900 nm & 930 nm SD-OCT Systems, Imperial	\$12,730.80	Lead Time

*1	Qty	Docs	Part Number - Universal		Price	Available
	<input type="text"/>		OCTG-900	Standard Scanner for 900 nm & 930 nm SD-OCT Systems	\$12,200.35	Lead Time
	<input type="text"/>		OCTH-900	Handheld Scanner for 900 nm & 930 nm SD-OCT Systems	\$10,502.91	Lead Time

*1	Qty	Docs	Part Number - Metric		Price	Available
	<input type="text"/>		OCTP-900/M	Customer Inspired! User-Customizable Scanner for 900 nm & 930 nm SD-OCT Systems, Metric	\$12,730.80	Lead Time

[Add To Cart](#)


Scan Lens Kit (Required OCT System Component)

- ▶ Telecentric Scan Lenses Provide a Flat Imaging Plane
- ▶ Lens AR Coated for 800 - 1100 nm
- ▶ Scan Lens Kits for Standard / User-Customizable Scanners Include
 - Telecentric Scan Lens
 - Illumination Tube
 - IR Card
 - Calibration Target
- ▶ Compact Scan Lens Kits Designed for the OCTH-900 Handheld Scanner with Integrated Scan Lens and Illumination Ring

**Also Required:
Base Unit
Scanning System**

To be functional, an OCT system build must include a base unit, a scanning system, and a scan lens kit.

The cross-section image below of a banana was taken with the OCT-LK3-BB scan lens kit using a Callisto Series OCT system. Choose a scan lens kit that provides the right resolution and focal length for your application.



Thorlabs' Scan Lens Kits enable easy exchange of scan lenses in an OCT system, providing the flexibility to tailor imaging resolution or working distance for each application. Based on our line of OCT telecentric scan lenses, these lens kits minimize image distortion without extensive post-image processing and maximize coupling of the light scattered or emitted from the sample surface into the detection system. As seen in the table below, we offer scan lens kits compatible with the standard (Item # OCTG-900) and user-customizable (Item # OCTP-900) scanners, as well as two lens kits compatible with the handheld scanner (Item # OCTH-900).

Each kit includes a telecentric scan lens, illumination tube, IR card, and calibration target. The included illumination tube serves as a light guide that channels light from the LED illumination ring down to the sample area. The IR card and calibration target are provided for calibration of the scanning mirror and lens kit, ensuring the best image quality when swapping between scan lenses.

Item #	OCT-LK2-BB	OCT-LK3-BB	OCT-LK4-BB	OCTH-LK20-BB	OCTH-LK30-BB
Click Image to Enlarge					
Design Wavelength	900 nm / 930 nm			900 nm / 930 nm	
Compatible Scanner	OCTG-900 (Standard) or OCTP-900 (User-Customizable)			OCTH-900 Handheld Scanner	
Lateral Resolution ^a	4 µm	8 µm	12 µm	9 µm	14 µm
Focal Length	18 mm	36 mm	54 mm	20 mm	30 mm
Working Distance	3.4 mm (with Tube) ^b 7.5 mm (without Tube)	24.9 mm (with Tube) ^b 25.1 mm (without Tube)	41.6 mm (with Tube) ^b 42.3 mm (without Tube)	12 mm	22 mm
Field of View	6 mm x 6 mm	10 mm x 10 mm	16 mm x 16 mm	Ø8 mm	Ø10 mm
Lens Threading	M25 x 0.75	M25 x 0.75	M25 x 0.75	M20 x 0.5 (For Z-Spacer)	M14 x 0.5 (For OCTH-900)

a. 1/e² Beam Diameter at Focus
b. The illumination tube is user-removable.

Based on your currency / country selection, your order will ship from Newton, New Jersey

*1	Qty	Docs	Part Number - Universal		Price	Available
	<input type="text"/>		OCT-LK2-BB	OCT Scan Lens Kit, 18 mm FL, 900 nm / 930 nm	\$2,056.02	Lead Time
	<input type="text"/>		OCT-LK3-BB	OCT Scan Lens Kit, 36 mm FL, 900 nm / 930 nm	\$1,406.75	Lead Time
	<input type="text"/>		OCT-LK4-BB	OCT Scan Lens Kit, 54 mm FL, 900 nm / 930 nm	\$1,406.75	Lead Time
	<input type="text"/>		OCTH-LK20-BB	OCT Scan Lens Kit for OCTH-900, 20 mm FL, 900 nm / 930 nm	\$1,071.29	Lead Time
	<input type="text"/>		OCTH-LK30-BB	OCT Scan Lens Kit for OCTH-900, 30 mm FL, 900 nm / 930 nm	\$1,071.29	Lead Time

[Add To Cart](#)

Reference Length Adapters (Required for Standard Scanners)

- ▶ Arm Adapters for Matching Reference Path Length to the Sample Path Length
- ▶ Use Multiple Reference Adapters for Rapid Switching Between Scan Lens Kits
- ▶ Must be Purchased with Standard Scanner (Item # OCTG-900)

Item # ^a	Compatible Scan Lens Kit
OCT-RA2	OCT-LK2-BB
OCT-RA3	OCT-LK3-BB



These adapters adjust the reference arm path length within the OCTG-900 Standard Scanner to match the sample path length of the scan lens used. Choose from three options that are compatible with the scan lens kits sold above. Reference length adapters also enable the user to quickly swap between different scan lens kits without going through extensive adjustments during each switch. The table to the right provides a compatibility list to help select the appropriate reference adapters.

OCT-RA4	OCT-LK4-BB
a. Multiple reference adapters can be purchased for rapid switching between scan lens kits.	

Based on your currency / country selection, your order will ship from Newton, New Jersey

+1	Qty	Docs	Part Number - Universal	Price	Available
<input type="text"/>	<input type="text"/>		OCT-RA2 Length Adapter for SD-OCT Standard Scanner & OCT-LK2(-BB) Scan Lens Kit	\$551.88	Lead Time
<input type="text"/>	<input type="text"/>		OCT-RA3 Length Adapter for SD-OCT Standard Scanner & OCT-LK3(-BB) Scan Lens Kit	\$551.88	Lead Time
<input type="text"/>	<input type="text"/>		OCT-RA4 Length Adapter for SD-OCT Standard Scanner & OCT-LK4(-BB) Scan Lens Kit	\$551.88	Lead Time

Sample Z-Spacers (Optional Accessories)



- ▶ Sample Z-Spacers Position Scanner at Optimal Working Distance From Sample
- ▶ Ring (Air) and Immersion (Liquid) Z-Spacers Available
- ▶ Two Z-Spacers Recommended for Use with OCTH-900 Handheld Scanner

Thorlabs offers both ring and immersion style sample Z-spacers that enable optimal positioning of a scanning system relative to the sample. The OCT-AIR3, OCT-IMM3, and OCT-IMM4 Z-Spacers feature knurled rings that allow the spacing distance to be adjusted and locked in place for increased stability. Several Z-spacer designs are available; please see the table below for compatibility with our scanners and lens kits.

Additionally, we offer two ring-style Z-spacers that are designed specifically for the OCTH-900 Handheld Scanner; these spacers greatly assist in maintaining the correct sample working distance when using the handheld scanner. The spacing distance on the OCTH-AIR20 and OCTH-AIR30 Z-Spacers can be adjusted by rotating the spacer.

[Click to Enlarge](#)
Z-Spacers for the OCTG-900 and OCTP-900(M) Scanners

Our ring-style Z-spacers provide a distance guide between the scanner and sample. The sample is in contact with the ring-shaped tip of the spacer and should only be used when air is the scanning medium. In contrast, our immersion spacers are equipped with a glass plate that contacts the sample surface within the scanning area. Unlike the ring-style spacers, immersion spacers enable access to samples contained within a liquid environment while also providing sample stabilization. Better index matching and a tilted glass plate also help reduce strong back-reflections from the sample surface and enhances the contrast of the image.



Item # ^a	Type	Adjustable	Adjustment Range	Lockable	Compatible Scanner	Compatible Scan Lens Kit
OCT-AIR3	Ring (Air)	Yes	+3.5 mm / -1.0 mm	Yes	OCTG-900 OCTP-900(M)	OCT-LK3-BB
OCT-IMM3	Immersion	Yes	-3.4 mm / -1.1 mm	Yes		OCT-LK4-BB
OCT-IMM4	Immersion	Yes	+1.0 mm / -17.0 mm	Yes	OCTH-900 ^a	OCTH-LK20-BB
OCTH-AIR20	Ring (Air)	Yes	±4 mm	No		OCTH-LK30-BB
OCTH-AIR30	Ring (Air)	Yes	±2 mm	No		OCTH-LK30-BB

a. We recommend purchasing a sample Z-spacer if using the OCTH-900 handheld scanner.

[Click to Enlarge](#)
Z-Spacers for the OCTH-900 Handheld Scanner

Based on your currency / country selection, your order will ship from Newton, New Jersey

+1	Qty	Docs	Part Number - Universal	Price	Available
<input type="text"/>	<input type="text"/>		OCT-AIR3 Ring-Style Sample Z-Spacer for OCT-LK3(-BB) Scan Lens Kit	\$772.63	Lead Time
<input type="text"/>	<input type="text"/>		OCT-IMM3 Immersion-Style Sample Z-Spacer for OCT-LK3(-BB) Scan Lens Kit	\$938.20	Lead Time
<input type="text"/>	<input type="text"/>		OCT-IMM4 Immersion-Style Sample Z-Spacer for OCT-LK4(-BB) Scan Lens Kit	\$1,048.57	Lead Time
<input type="text"/>	<input type="text"/>		OCTH-AIR20 Ring-Style Sample Z-Spacer for OCTH-LK20(-BB) Scan Lens Kit	\$199.11	Lead Time
<input type="text"/>	<input type="text"/>		OCTH-AIR30 Ring-Style Sample Z-Spacer for OCTH-LK30(-BB) Scan Lens Kit	\$199.11	Lead Time


Scanner Stand (Optional Accessory)



- ▶ Recommended Stand for Mounting Standard or User-Customizable Scanners
- ▶ Focus Block with Coarse/Fine Z-Axis Travel on Ø1.5" Stainless Steel Post
- ▶ 12" x 14" (300 mm x 350 mm) Aluminum Breadboard with 1/4"-20 (M6) Tapped Holes

For convenient mounting of our Standard or User-Customizable Scanners, we offer a scanner stand that is ideal for use in vibration-sensitive studies such as angiography. It consists of a post-mounted focus block with knobs that provide both coarse (40 mm/rev) and fine (225 µm/rev) z-axis travel. A rotation and height collar underneath the focus block allows it to rotate 45° in order to move the scanner head away from the sample to make adjustments.

[Click to Enlarge](#)
Translation Stage (sold below) to the OCT-STAND and the OCT-XYR1/M Translation Stage to the OCT-STAND/M directly underneath the scan lens. A 1/4"-20 (M6) counterbore is also provided for securing the Ø1.5" post.



[Click for Details](#)
The focus block can be rotated 45° to move the scanner head away from the sample.

Based on your currency / country selection, your order will ship from Newton, New Jersey

+1	Qty	Docs	Part Number - Imperial	Price	Available
<input type="text"/>	<input type="text"/>		OCT-STAND Stand for Standard and User-Customizable OCT Scanning Systems, 1/4"-20 Tapped Holes	\$2,164.24	Today

+1	Qty	Docs	Part Number - Metric	Price	Available
<input type="text"/>	<input type="text"/>		OCT-STAND/M Stand for Standard and User-Customizable OCT Scanning Systems, M6 Tapped Holes	\$2,164.24	Today

Translation Stage (Optional Accessory)



- ▶ Optional Translation Stage with 0.5" (13 mm) of XY Travel and 360° Rotation
- ▶ Includes Cover Plate for Sample Mounting
- ▶ Can Mount Optomechanics by Removing Cover Plate

[Click to Enlarge](#)

Precise translation and rotation are often required for optimal positioning of a sample before and during OCT imaging. The OCT-XYR1(M) is an XY linear translation stage with a rotating platform and solid plate for sample mounting and easy cleaning. The OCT-XYR2 or OCT-XYR1/M stage can be secured to the OCT-STAND or OCT-STAND/M, respectively, using the 1/4" (M6) counterbores at the corners. The top plate is removable for access to 4-40, 8-32 (M4) and 1/4"-20 (M6) tapped holes and an SM1-threaded (0.025"-40) central hole for mounting optomechanical components. The XYR1 Solid Sample Plate can be purchased separately as a direct replacement for the top plate.

Specifications	
Horizontal Load Capacity (Max)	10 lbs (4.5 kg)
Mounting Platform Dimensions	Ø4.18" (Ø106 mm)
Stage Height	1.65" (41.8 mm)
Linear Translation Range	1/2" (13 mm)
Travel per Revolution	0.025" (0.5 mm)
Graduation	0.001" (10 µm) per Division



[Click to Enlarge](#)
The cover plate is removable for access to tapped holes and the SM1-threaded central hole.

The X and Y micrometers offer 1/2" (13 mm) of travel with graduations every 0.001" (10 µm). The stage's rotation and translation can be freely changed without compromising the stability of attached components. An engraved angular scale along the outer edge of the stage's rotating platform allows the user to set the angular orientation of the stage, which can then be fixed using the 5/64" (2 mm) hex locking set screws. Locking the rotation of the stage does not prevent XY translation using the actuators.

Based on your currency / country selection, your order will ship from Newton, New Jersey

+1	Qty	Docs	Part Number - Imperial	Price	Available
<input type="text"/>	<input type="text"/>		OCT-XYR1 XY Stage with Solid Top Plate, 1/2" Travel, 360° Rotation, Imperial Taps	\$772.63	Today

+1	Qty	Docs	Part Number - Metric	Price	Available
<input type="text"/>	<input type="text"/>		OCT-XYR1/M XY Stage with Solid Top Plate, 13 mm Travel, 360° Rotation, Metric Taps	\$772.63	5-9 Days

Additional OCT Imaging Systems & Components

OCT Systems Comparison	Telesio Series OCT Systems and Components	Publications
OCT Systems for OEM Applications	Vega Series OCT Systems and Components	Application Articles
Scanmedx Series OCT Systems and Components	Telesio Series PR-OCT Systems and Components	
Callisto Series OCT Systems and Components	OCT Components	