

TRUEplus™ Urine Test Strips for Ketones

Instructions for Use

Summary: In the body, carbohydrates are converted to glucose.

Glucose is the body's primary source of energy. Insulin is needed to help process glucose in blood to supply the body with energy. When blood glucose is not available for energy, the body will use fat stores to produce ketones for energy. Excess ketones are discarded in urine.¹

For low carbohydrate dieters who are not diabetic, low intake of complex carbohydrates and sugars promotes use of ketones from fat stores rather than blood glucose as the primary source of energy for the body. When this occurs, the body produces a steady state of ketones (ketosis).¹

For diabetics, a lack of sufficient insulin prevents the body from using blood glucose properly. Without enough glucose for energy, the body produces ketones from fat and muscle for energy (ketoacidosis),

Expected
negative (0)
gives negative
Healthcare
Large (160 m

Unusual Tes

1. Check 'Use
or it is 2 m
with Strips
2. Check for di
if vial cap wa
heat or cold.

If you have que

Storing Test S

- Test Strips

Contact with

TRUEplus

KETONE

test strips

(Acetoacetic Acid)

Reagent Strips

for Urinalysis

- 15 seconds test time
- Convenient flip-top vial
- 10

Store in a dry place at room

temperature below 86°F (30°C).

Manufactured by



TRUEplus

KETONE

test strips

Acetoacetic Acid
Reagent Strips
for Urinalysis

Store in a dry place at room
temperature below 86°F (30°C).
DO NOT REFRIGERATE OR FREEZE.

Manufactured

RUEPLUS™

ETONE

Fast strips

dry place at room

below 86°F (30°C).

GENERATE OR FREEZE.

ay from sunlight.