
TRACERS®

Bone Cement Opacifier
Opacifiant de ciment osseux
Opacizzante del cemento osseo
Opakmittel für Knochenzement
Opacador de cemento óseo
Opacificador de cimento ósseo
Troebelingsmiddel voor botcement
Knoglecementopacifikator
Luusementin samennin
Opakmedel för bencement
Μέσο αδιαφανοποίησης οστικού τσιμέντου

Parallax®

English Instructions for Use

INDICATIONS FOR USE:

Parallax TRACERS Bone Cement Opacifier is intended for use as an additive to Secour Acrylic Resin to provide radiopacity for imaging purposes.

CAUTION:

Federal (USA) law restricts this device to sale by or on the order of a physician. Radiopacity may mask underlying pathological conditions. Radiopacity may also make it difficult to assess radiographically the ingrowth of new bone.

CONTRAINDICATIONS:

Adhere to contraindications provided in the Secour Instructions for Use.

CONTENTS:

Each vial contains 5g Barium Sulfate.

PROCEDURE:

Preparation of the Mix For Repairing Cranial Defects.

1. Remove the TRACERS vial from the pouch using standard sterile technique.
2. Empty the Secour powder into the TRACERS vial.
3. Seal the vial and shake to ensure the powder is homogenous.
4. Continue mixing according to the Secour Instructions for Use.
5. When the mixture has reached dough stage, mix again thoroughly to disperse the particles evenly in the dough mass.

Note: Addition of TRACERS particles may increase Secour dough time by up to 1 minute, may increase Secour working time by up to 1 minute, and may increase the setting time by up to 3 1/2 minutes.

MANIPULATION AND PLACEMENT:

1. Apply the mix to the defect site per Secour Instructions for Use.

STERILIZATION:

This product is sterile and non-pyrogenic if packaging is not opened, damaged, or broken. Single use only.

HOW SUPPLIED:

The device is provided sterile and intended for **SINGLE USE ONLY. DO NOT clean, resterilize, or reuse, as this may result in product malfunction, failure, or patient injury, and which may also expose the patient to transmitted infectious disease.** Refer to package labeling for sterilization method used.

PRODUCT INFORMATION DISCLAIMER:

ArthroCare Corporation, excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to any implied warranties of merchantability or fitness, since handling and storage of this device by the user, as well as factors relating to the patient's diagnosis, treatment and other matters beyond ArthroCare Corporation's, direct control affect this device and the results obtained with the use. ArthroCare Corporation, shall not be liable for any incidental, or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. ArthroCare Corporation, neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

RETURN CONDITIONS:

In the event the device must be returned for any reason, return the unopened material in its original package and shipping box. Call +1-800-797-6520 to receive a return goods authorization number prior to return shipment.

© 2004 ArthroCare Corporation, ArthroCare, Parallax, TRACERS, Clear-View, Core-Assure, EZ Flow, Secour and CDS are trademarks and/or registered trademarks of ArthroCare Corporation.

US patents 6019776, 6033411, 6231615, 6309420, 6348055, 6383190 and patents pending.

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