




CLIA Waived



iCup[®] Dx Drug Screen

Rapid Drug Screening Made Simple.


- CLIA Waived
- Up to 12 drugs
- Simple Process - tilt, read at 5 minutes
- All Inclusive - no urine exposure
- Leak-proof, secure collection cup
- Tester activates the test when ready
- Integrated Temperature Strip
- Tamper Resistant



Innovative Drug Screening Solutions
Quality • Choices • Support

Instant Technologies, Inc.




an  inverness medical company

iCup[®] Dx Drug Screen

Rapid Drug Screening Made Simple.



Detect the presence of drugs and adulterants in minutes.
It's drug screening efficiency you can count on.




iCup Dx 12-Panel with SVT


Drug Class:	Abbreviation:	Cutoff Levels:
Amphetamine.....	AMP.....	1000ng/ml
Barbiturates.....	BAR.....	300ng/ml
Benzodiazepines.....	BZO.....	300ng/ml
Cocaine.....	COC.....	300ng/ml
Ecstasy.....	MDMA.....	500ng/ml
Marijuana.....	THC.....	50ng/ml
Methadone.....	MTD.....	300ng/ml
Methamphetamine.....	mAMP.....	1000ng/ml
Opiates.....	OPI.....	2000ng/ml
Oxycodone.....	OXY.....	100ng/ml
Phencyclidine.....	PCP.....	25ng/ml
Tricyclic Antidepressants.....	TCA.....	1000ng/ml

Instructions

- 1 Collect specimen and secure test lid**
- 2 Tilt cup on its side to activate test**
- 3 Read the results**



iCup Dx without SVT.



iCup Dx with SVT.



Read results at five (5) minutes.



NEGATIVE NON-NEGATIVE INVALID

	ADULTERATION COLOR COMPARISON CHART		
	ABNORMAL (LOW)	NORMAL	ABNORMAL (HIGH)
Creatinine (CR) 45 seconds	<10 10	20 50 100 200 mg/dl	≥150 mg/dl
Nitrite (Ni) 45 seconds		0 0.1-0.2 0.5-5.0	≥10
pH Immediate	2 3	4 5 7 8	
Bleach (BI) 45 seconds		Negative	Positive
Specific Gravity (S.G.) 80 seconds	1.000	1.005 1.015 1.025	≥1.030

SAMPLE ONLY - Color chart NOT Exact

Catalog #

Drugs/Panel

Optional Specimen Validity Tests (SVT)

I-DXA-187-016 *	8	COC, THC, OPI, AMP, mAMP, BZO, OXY, MTD (Cr, Ni, pH, BI, S.G.)
I-DXA-1107-015 *	10	COC, THC, OPI, AMP, MDMA, PCP, BAR, BZO, MTD, TCA (Cr, Ni, pH, BI, S.G.)
I-DXA-1107-142 *	10	COC, THC, OPI, AMP, mAMP, MDMA, BZO, BAR, OXY, MTD (Cr, Ni, pH, BI, S.G.)
I-DXA-1127-023	12	COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, OXY, MTD, TCA, MDMA (Cr, Ni, pH, BI, S.G.)

* Coming Soon



*A rapid, all-inclusive drug screen with optional
Specimen Validity Tests.*

(BZO/BAR)+(COC/THC)+(MET1000/OPI2000)
+(MTD/TCA)+(OXY/MDMA)+(PCP/AMP)

Item# I-DXA-1127-023

Lot: 90856 **Exp:** 12/2010

Contents: Individually packaged devices (25), package insert (1)
Security Seals (25), Color Comparison Chart (1)

For in vitro diagnostic use only
Storage temperature: 2 - 30°C (36 - 86°F)

Manufactured for:  **inverness medical**
Orlando, FL 32810 USA

made in the USA

iCup[®]

(mAMP/COC/THC)+(AMP/OPI)+
(OXY/PPX)+(BZO/BAR)+<MDMA>

Catalog Number: i-DOA-1107-051

Lot: DOA3030361 Exp: 2014/12

Manufactured for Instant Technologies, Inc., Norfolk, VA23502

For *in vitro* diagnostic use only
Store at any temperature from 2 - 30°C (36 - 86°F)

L020453-02

iCup[®]

25
CUPS

A rapid, one step screening test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine.

**(mAMP/COC/THC)+(AMP/OPI)+
(OXY/PPX)+(BZO/BAR)+<MDMA>
Catalog Number: i-DOA-1107-051**

Lot:DOA3030361

Exp:2014/12

For in vitro diagnostic use only

Store at any temperature from 2 - 30°C (36 - 86°F)

This assay provides only a preliminary result. A quantitative analytical method is needed to obtain a confirmed result (see package insert).

Manufactured for Instant Technologies, Inc., Norfolk, VA 23502

L040437-01



Instant Technologies One Step Multi-Drug Screen Test Card with the Integrated *iCup® iCup®*

Instruction Sheet for testing of any combination of the following drugs:
AMP/BAR/BUP/BZO/COC/THC/MTD/MDA/MAMP/IND/MA/MP/OP/OX/XY/PC/P/PP/PT/CA

Available with Specimen Validity Tests (S.V.T.) for Oxidant/POCC, Specific Gravity, pH, Nitrite, Glutaraldehyde and Creatinine

A rapid, one step screening test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine.
For healthcare professionals including professionals at point of care sites.
Immunoassay for *in vitro* diagnostic use only.

INTENDED USE

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} is a lateral flow immunographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in urine at the following cut-off concentrations:

Drug	Calibrator	Cut-off
Amphetamine (AMP 1,000)	4-N-Acetylmethamphetamine	1,000 ng/mL
Amphetamine (AMP 300)	4-Amphetamine	300 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	10 ng/mL
Buprenorphine (BUP)	Buprenorphine	300 ng/mL
Cocaine (COC 300)	Benzoylcegonine	150 ng/mL
Cocaine (COC 150)	11-nor-2'-THC-COOH	30 ng/mL
Marijuana (MTD)	11-nor-2'-THC-COOH	300 ng/mL
Methamphetamine (MAMP 1,000)	4-Methamphetamine	1,000 ng/mL
Methamphetamine (MAMP 500)	Methamphetamine (MAMP 500)	500 ng/mL
Methylenedioxymethamphetamine (MDMA)	4-Methylenedioxymethamphetamine	500 ng/mL
Opiate (MOP 300)	Morphine	300 ng/mL
Opiate (MOP 150)	Morphine	150 ng/mL
Oxycodone (OCX)	Oxycodone	100 ng/mL
Phenothiazine (PPH)	Phenothiazine	25 ng/mL
Propoxyphene (PPX)	Propoxyphene	300 ng/mL
Tetrahydrocannabinol (THC)	THC (delta-9-tetrahydrocannabinol)	1,000 ng/mL

Configurations of the One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} come with any combination of the drugs listed above. This assay provides only a preliminary analytical test result. A more specific alternate chemical method should then be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred method for final identification and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

SUMMARY

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs in urine.

AMPHETAMINE (AMP 1,000)

Amphetamine is a Schedule II controlled substance available by prescription (Dexedrine®) and is also available on the illicit market. Amphetamine is a class of potent sympathomimetic agents with therapeutic applications. They are chemically related to the human body's natural catecholamines, epinephrine and norepinephrine. Acute higher doses lead to enhanced stimulation of the central nervous system (CNS) and induce euphoria, alertness, increased energy and power. Cardiovascular responses to amphetamines include increased blood pressure, increased heart rate, increased pulse rate, decreased appetite, anorexia, hallucinations, and psychotic behavior. The effects of amphetamines are acute responses that occur in the urine and the drug has a half-life of 4-24 hours in the body. About 30% of amphetamines are excreted in the urine in unchanged form, with the remainder as hydroxylated and deaminated derivatives.

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of amphetamines in urine exceeds 1,000 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).

AMPHETAMINE (AMP 300)

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when amphetamines in urine exceed 300 ng/mL. See AMPHETAMINE (AMP 1,000) for the summary.

BARBITURATES (BAR)

Barbiturates are CNS depressants. They are used therapeutically as sedatives, hypnotics, and anesthetic agents. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of barbiturates leads to tolerance and physical dependence. Short-acting barbiturates taken at 400 mg/day for 2-3 months can produce a clinically significant degree of physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death.

Only a small amount (less than 5%) of most barbiturates are excreted unaltered in the urine. The approximately 40% of barbiturates are excreted unaltered in the urine.

Short acting (e.g. Secobarbital) 100 mg PO (oral)
Long acting (e.g. Phenobarbital) 400 mg PO (oral)
4.5 days
7 days

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of barbiturates in urine exceeds 300 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for Barbiturate positive specimens.

BENZODIAZEPINES (BZO)

Benzodiazepines are medications that are frequently prescribed for the symptomatic treatment of anxiety and insomnia. They induce their effects via specific receptors involving a neurochemical called gamma aminobutyric acid (GABA). Benzodiazepines have been found to be highly effective, but have been replaced in some surgical and medical procedures, and for the treatment of benzodiazepines are also used as sedatives before some surgical and medical procedures, and for the treatment of benzodiazepines are also used as sedatives before some surgical and medical procedures. Risk of physical dependence increases if benzodiazepines are taken regularly (e.g. daily) for a few months, especially at higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal upset, feeling unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception.

Only trace amounts (less than 1%) of most benzodiazepines are excreted unaltered in the urine; most of the concentration in urine is conjugated drug. The detection period for benzodiazepines in urine is 3-7 days. The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of benzodiazepines in urine exceeds 300 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for benzodiazepine positive specimens.

BUPRENORPHINE (BUP)

Buprenorphine is a semi-synthetic opioid analgesic often used in the treatment of opioid addiction. The drug is sold under the trade name Subutex®. Buprenorphine is a partial agonist at the mu-opioid receptor and is chemically related to morphine, but with a much longer duration of action. Buprenorphine, in combination with Naloxone HCl Therapeutically Buprenorphine, is a form of medical care offered to opioid addicts. Substitution treatment is a form of medical care offered to opioid addicts (formerly, Buprenorphine is as effective as Methadone but demonstrates a lower level of physical dependence. Concentrations of free Buprenorphine and Nalbupropine in urine may be less than 1 ng/ml after therapeutic administration, but can range up to 20 ng/ml in abuse situations. The plasma half life of Buprenorphine is 2-4 hours. While complete elimination of a single dose of the drug can take as long as 6 days, the window of detection for the parent drug in Substituted Abuse of Buprenorphine is approximately 9 days.

Buprenorphine has been reported in many countries where various forms of the drug are available. The drug has been abused from intranasal, intranasal, intranasal and intravenous (IV) routes. Buprenorphine has been abused via intravenous, sublingual, intranasal and intravenous routes. The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of Buprenorphine in urine exceeds 10 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for buprenorphine positive specimens.

COCAINE (COC 300)

Cocaine is a potent central nervous system stimulant and a local anesthetic. Initially, it brings about extreme energy and restlessness while gradually resulting in tremors, over-sensitization and spasms. In large amounts, cocaine causes fever, unresponsiveness, difficulty in breathing and uncoordinatedness.

Cocaine is often self-administered by nasal inhalation, intravenous injection and free-base smoking. It is excreted in the urine in a short time primarily as benzoylcegonine. Benzoylcegonine, a major metabolite of cocaine, has a half-life of 3-5 hours and can generally be detected for 24-48 hours after cocaine exposure (6-8 hours) than cocaine (0.5-1.5 hours), and can generally be detected for 2-48 hours after cocaine exposure (6-8 hours) than cocaine (0.5-1.5 hours), and can generally be detected for 2-48 hours after cocaine exposure (6-8 hours) than cocaine (0.5-1.5 hours), and can generally be detected for 2-48 hours after cocaine exposure (6-8 hours) than cocaine (0.5-1.5 hours).

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of benzoylcegonine in urine exceeds 300 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).

COCAINE (COC 150)

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of benzoylcegonine in urine exceeds 150 ng/mL. See COCAINE (COC 300) for the summary.

MARIJUANA (THC)

THC (delta-9-tetrahydrocannabinol) is the primary active ingredient in cannabis (marijuana). When smoked or orally administered, THC produces euphoric effects. Users have impaired short-term memory, altered perceptions, and they may also experience transient episodes of confusion and anxiety. Long-term, reliably heavy use may be associated with behavioral disorders. The peak effect of marijuana administered by smoking occurs in 20-30 minutes and the duration is 90-120 minutes after one cigarette. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 3-10 days after smoking. The main metabolite excreted in the urine is 11-nor-9-tetrahydrocannabinol-9-carboxylic acid (THC-COOH).

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of THC-COOH in urine exceeds 300 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).

METHADONE (MTD)

Methadone is a narcotic analgesic prescribed for the management of moderate to severe pain and for the treatment of opiate dependence (Nexin, Viodin, Percopin, morphine). The pharmacology of oral methadone is very different from methadone. Oral methadone is partially passed in the liver for later use. IV methadone acts more like heroin. In most patients, methadone is long acting pain reliever, producing a long duration of action. When given orally, methadone treats the acute pain of withdrawal. Methadone is a long acting pain reliever that can be prescribed to patients with the emotional roller coaster that most opiates produce. Methadone, if taken for long periods and at large doses, can lead to a very long withdrawal period. The withdrawals from methadone are more prolonged and troublesome than those provoked by heroin cessation, yet the substitution and phased removal of methadone is an acceptable method of detoxification for patients and therapists.

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of methadone in urine exceeds 300 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for methadone positive specimens.

METHAMPHETAMINE (MAMP 1,000)

Methamphetamine is an addictive stimulant drug that strongly activates certain systems in the brain. Methamphetamine is closely related chemically to amphetamine, but the CNS effects of methamphetamine are greater. Methamphetamine is used by drug abusers and has a high potential for abuse and dependence. The drug can be taken orally, injected, or smoked. Methamphetamine is used to enhance stimulation of the CNS and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to methamphetamine include increased blood pressure and cardiac arrhythmias. Methamphetamine responses produce anxiety, paranoia, hallucinations, psychotic behavior, and eventually, depression and exhaustion.

The effects of methamphetamine generally last 2-4 hours and the drug has a half-life of 9-24 hours in the body. Methamphetamine is excreted in the urine as amphetamine and oxidized and deaminated derivatives. However, 10-30% of methamphetamine is excreted unchanged. Thus, the presence of the parent compound in the urine indicates methamphetamine use. Methamphetamine is generally detectable in the urine for 3-5 days, depending on urine pH level. The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of methamphetamine in urine exceeds 1,000 ng/mL.

METHAMPHETAMINE (MAMP 500)

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of methamphetamine in urine exceeds 500 ng/mL. See METHAMPHETAMINE (MAMP 1,000) for the summary.

METHYLENEDIOXYMETHAMPHETAMINE (MDMA)

Methylenedioxymethamphetamine (ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity. Those who take the drug frequently report adverse effects, including increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release 5-HT, although the general opinion is that this is a secondary effect of the drug (Nichols and Oberlander, 1990). The one major effect of MDMA, occurring in virtually all people who took a reasonable dose of the drug, was to produce a subjective effect of MDMA, occurring in virtually all people who took a reasonable dose of the drug.

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of Methylenedioxymethamphetamine in urine exceeds 500 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for Methylenedioxymethamphetamine positive specimens.

OPIATE (MOP 300)

Opiate refers to any drug that is derived from the opium poppy, including the natural products, morphine and codeine, and the semi-synthetic drugs such as heroin. Opioid is more general, referring to any drug that acts on the opioid receptor.

Opioid analgesics comprise a large group of substances which control pain by depressing the CNS. Large doses of morphine can produce high tolerance levels, and also the dependency in users, and may lead to substance abuse. Morphine is excreted unmetabolized, and is also the major metabolic product of codeine and heroin.

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of morphine in urine exceeds 300 ng/mL.

OPIATE (OPI 2,000)

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of morphine in urine exceeds 2,000 ng/mL. See OPIATE (MOP 300) for the summary.

OXYCODONE (OXY)

Oxycodone is a semi-synthetic opioid with a structural similarity to codeine. The drug is manufactured by modifying thebaine, an alkaloid found in the opium poppy. Oxycodone, like all opiates agonists, provides pain relief by acting on opioid receptors in the spinal cord, brain, and possibly directly in the affected tissues. Oxycodone is prescribed for the relief of moderate to high pain under the well-known pharmaceutical trade names of OxyContin®, Tylox®, Percocet® and Percocet®. While Tylox, Percocet and Percocet contain only oral doses of oxycodone hydrochloride combined with other analgesics such as acetaminophen or aspirin, OxyContin consists of oxycodone hydrochloride in a time-release form. Oxycodone is known hydrophilic and is excreted in the urine as oxycodone and nortoxycodone. In a 24-hour urine, 33-61% of a single, 5 mg oral dose is excreted with 10-15% of the dose being unchanged drug (13-15%). The urine is expected to be similar to that of other opiates such as morphine.

The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of oxycodone in urine exceeds 100 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for oxycodone positive specimens.

PHENCYCLIDINE (PCP)

Phencyclidine, also known as PCP or Angel Dust, is a hallucinogen that was first marketed as a surgical anesthetic in the 1950's. It was removed from the market because patients receiving it became delirious and experienced hallucinations.

PCP is used in powder, capsule, and tablet form. The powder is either snorted or smoked after mixing it with a substance like sugar. PCP is most commonly administered by inhalation but can be used intravenously, intranasally, or rectally. PCP acts as a central nervous system stimulant and acts wily and unpredictable mood swings from euphoria to depression. Specific behavioral responses of PCP use include: increased heart rate and blood pressure. PCP can be found in urine within 4 to 6 hours after use and will remain in urine for 7 to 14 days, depending on factors such as metabolic rate, user's age, weight, activity, and diet. PCP is excreted in the urine as an unchanged drug (4% to 19%) and conjugated metabolites (25% to 30%). The One Step Multi-Drug Screen Test Card with the Integrated iCup® iCup®^{AB} yields a positive result when the concentration of phencyclidine in urine exceeds 25 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).

PROPOXYPHENE (PPX)

Propoxyphene (PPX) is a narcotic analgesic compound bearing structural similarity to methadone. As an analgesic, propoxyphene can be from 50-75% as potent as oral codeine. Darvocet™, one of the most common brand names for the drug, contains 50 mg of propoxyphene napsylate and 325-650 mg of acetaminophen. Propoxyphene is a narcotic analgesic with a half-life of approximately 1 to 2 hours post dose. In the case of overdose, propoxyphene blood concentrations can reach significantly higher than therapeutic levels. Propoxyphene has a longer half-life (30 to 36 hours) than parent propoxyphene (8 to 12 hours). The accumulation of nortropoxyphene has been with repeated doses may be largely responsible for resultant toxicity. The One Step Multi-Drug Screen Test Card with the Integrated iCup[®]/iCup[®]AD yields a positive result when the test strip is immersed in a urine specimen that exceeds 300 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for propoxyphene positive specimens.

TRICYCLIC ANTIDEPRESSANTS (TCA)

TCA (Tricyclic Antidepressants) are commonly used for the treatment of depressive disorders. TCA overdoses can result in profound CNS depression, cardiotoxicity and anticholinergic effects. TCA overdose is the most common cause of death from prescription drugs. TCAs are taken orally or sometimes by injection. TCAs are metabolized in the liver. Both TCAs and their metabolites are excreted in urine mostly in the form of metabolites for up to ten days. The One Step Multi-Drug Screen Test Card with the Integrated iCup[®]/iCup[®]AD yields a positive result when the test strip is immersed in a urine specimen that exceeds 100 ng/mL. At present, the Substance Abuse and Mental Health Services Administration (SAMHSA) does not have a recommended screening cut-off for tricyclic antidepressant positive specimens.

S.V.T. SUMMARY

(Information regarding Specimen Validity Tests does not require FDA review.) The strip contains chemically treated reagent pads, 3-5 minutes following the activation of the reagent pads by the addition of a urine specimen. The color change in the test strip is compared with the printed color chart card. The color comparison provides a semi-quantitative result. The test is used to determine certain urinary characteristics such as pH and specific gravity, pH, nitrite, glutaraldehyde, and creatinine in human urine which can help assess the integrity of the urine sample.

WHAT IS ADULTERATION?

Adulteration is the tampering of a urine specimen with the intention of altering the test results. The use of adulterants present in the urine, dilution may also be used to alter the results. Adulterants may be used to alter the results of the test. One of the best ways to test for adulteration or dilution is to determine certain urinary characteristics such as pH and specific gravity and to detect the presence of oxidants/PCC, specific gravity, pH, nitrite, glutaraldehyde and creatinine in urine.

- Oxidants/PCC (Pyridinium chlorochromate) tests for the presence of oxidizing agents such as bleach and adulterants. Normal human urine should not contain oxidants or PCC.
- Specific gravity tests for sample dilution. The normal range is from 1.003 to 1.030. Values outside this range may be the result of specimen dilution or adulteration.
- pH tests for the presence of acidic or alkaline adulterants in urine. Normal pH levels should be in the range of 5.0 to 8.0. Values outside of this range may indicate the sample has been altered.
- Nitrite tests for the presence of nitrite adulterants in urine. Normal urine should contain no nitrite. Preliminary positive results generally indicate the presence of an adulterant.
- Glutaraldehyde tests for the presence of an aldehyde. Adulterants such as UrinAM and Clear Choice used in some urine specimen vials may cause false negative screening results by disrupting the enzyme used in certain urine specimen vials. Glutaraldehyde is not normally found in urine; therefore, detection of glutaraldehyde in a urine specimen is a waste product of creatine, an amino-acid contained in muscle tissue and found in urine. A person may attempt to foil a test by drinking excessive amounts of water or diuretics such as herbal teas to "flush" the system. Creatinine and specific gravity are two ways to check for dilution and flushing, which are the most common mechanisms used in an attempt to circumvent drug testing. Low creatinine and specific gravity may indicate dilute urine. The absence of creatinine (< 5 mg/dl) is indicative of a specimen not consistent with human urine.

PRINCIPLE

The One Step Multi-Drug Screen Test Card with the Integrated iCup[®]/iCup[®]AD is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen, will compete with the drug-protein conjugate and a visible colored line will show up in the test region. The antibody will then react with the drug above the cut-off concentration will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test region.

A drug-positive urine specimen will not generate a colored line in the specific test region of the strip because of drug competition, while a drug-negative urine specimen will generate a line in the test region because of the test strip. To serve as a procedural control, a colored line will always appear at the control region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

Each test contains anti-drug mouse monoclonal antibody and corresponding drug-protein conjugates. The control line contains goat anti-rabbit IgG polydonal antibodies and rabbit IgG.

S.V.T. REAGENTS

Adulteration Pad	Reactive Indicator	Buffer and Inhibitory Ingredients
Oxidants/PCC	0.25%	99.75%
Specific Gravity	0.25%	99.75%
pH	0.05%	99.95%
Nitrite	0.07%	99.93%
Glutaraldehyde	0.02%	99.98%
Creatinine	0.04%	99.96%

PRECAUTIONS

- For healthcare professionals including professionals at point of care sites.
- Immunoassay for in vitro diagnostic use only. Do not use after the expiration date.
- All specimens should be considered potentially infectious.
- The used test cup should be discarded according to federal, state and local regulations.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C (36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test devices must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

The urine specimen must be collected in a clean, dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should not be thawed and retested within 24 hours. When tests include S.V.T., storage of urine specimens should not exceed 2 hours at room temperature or 4 hours refrigerated prior to testing. For best results, test specimens immediately following collection.

MATERIALS Provided

- Cups with multi-drug panels [Note: A Fahrenheit temperature strip is affixed to aid in the determination of testing validity. Please use this temperature strip in conjunction with your Drug Free Policy (if applicable)].
- Adhesive Seal Label (if applicable)
- Security seal label
- Package insert
- Procedure card


Materials Required But Not Provided

- A timer or any kind of timing device such as a wrist watch is required to run the test.
- External controls


DIRECTIONS FOR USE

- Allow the test cup, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.
- Bring the pouch to room temperature before opening it. Remove the cup from the sealed pouch and use it as directed.
- Donor provides specimen.
- Technician replaces and secures cap while the cup is on a flat surface.
- Donor dates and initials the security seal and attaches the security seal over the cup cap.
- Technician peels off label to reveal adulteration strip(s), if applicable.
- Technician peels off the label on the multi-drug test card to view results.
- The adulteration strip(s), if applicable, should be read between 3-5 minutes. Compare the colors on the adulteration strip to the color chart. If the results indicate adulteration, do not read the drug test results. If adulteration is not indicated, the drug test results are stable for up to sixty minutes. For detailed operation instructions, please refer to the Procedure Card and Color Chart.
- If preliminary positive results are observed, please send the cup to the laboratory for confirmation.

Cup without SVT



Cup with SVT



INTERPRETATION OF RESULTS

NEGATIVE: A colored line appears in the Control region (C) and a colored line appears in the Test region (Drug/T) next to a specific drug tested. This negative result means that the drug concentrations in the urine sample are below the designated cut-off levels for a particular drug tested.

NOTE: The shade of the colored lines in the Test region may vary. The result should be considered negative whenever there is even a faint colored line.

POSITIVE: A colored line appears in the Control region (C) and NO line appears in the Test region (Drug/T) next to the name of a specific drug tested. The positive result means that the drug concentration in the urine sample is greater than the designated cut-off for a specific drug.

INVALID: NO line appears in the Control region (C), insufficient specimen volume or incorrect procedural techniques are used. If the result is still invalid, contact your distributor. Repeat the test with a new test cup. If the result is still invalid, contact your distributor.

SVT/ADULTERANT INTERPRETATION

Please refer to the color chart. Semi-quantitative results are obtained by comparing the reacted color blocks on the adulteration strips to the printed color blocks on the color chart. No color change is observed.

QUALITY CONTROL

A procedural control is included in the test. A line appearing in the Control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The One Step Multi-Drug Screen Test Card with the Integrated iCup[®]/iCup[®]AD provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{1,2,3}
- There is a possibility that technical or procedural errors, as well as interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result does not indicate level or intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- This test does not distinguish between drugs of abuse and certain medications.
- A positive test result may be obtained from certain foods or food supplements.

S.V.T.-ADULTERATION LIMITATIONS

- The adulteration tests, included with this product, are meant to aid in the determination of abnormal specimens. While comprehensive, these tests are not meant to be an alcoholic representation of possible adulterants.
- Oxidants/PCC: Normal human urine should not contain oxidants or PCC. The presence of high levels of oxidants in the specimen, such as ascorbic acid, may result in false negative results for the oxidants/PCC pad.
- Specific Gravity: Elevated levels of protein in urine may cause abnormally high specific gravity values.
- Nitrite: Nitrite is not a normal component of human urine. However, nitrite found in urine may indicate urinary tract infections or bacterial infections. Nitrite levels of > 20 mg/dl, may produce false preliminary positive glutaraldehyde results.
- Glutaraldehyde: Is not normally found in urine. However certain metabolic abnormalities such as ketoadcrosis (fasting, uncontrolled diabetes or high-protein diets) may interfere with the test results.
- Creatinine: Normal creatinine levels are between 20 and 350 mg/dL. Under rare conditions, certain kidney diseases may show dilute urine.

PERFORMANCE CHARACTERISTICS

ACCURACY

A side-by-side comparison was conducted using the One Step Multi-Drug Screen Test Card with the Integrated iCup[®]/iCup[®]AD and commercially available drug rapid tests. Testing was performed on 300 specimens per drug type previously collected from subjects present for drug screen testing. Presumptive positive results were confirmed by GC/MS. The following compounds were queried by GC/MS and contributed to the total amount of drugs found in presumptive positive urine samples tested.

Test	Compounds Contributing to GC/MS Totals
AMP	Amphetamine
BAR	Secobarbital, Barabital, Phenobarbital, Pentobarbital
BUP	Buprenorphine
EZO	Oxazepam, Nortriptyline, o-Hydroxyflazepam, Desmethylflazepam
COC	Benzoylcocaine
THC	11-nor- Δ^9 -tetrahydrocannabinol-9-carboxylic acid
MTPD	Mepharmidone
mAMP	Mepharmidone
MDMA	d-Methylenedioxymethylamphetamine
COB	Norphine, Codeine
COX	Oxycodone
PRC	Propylthiouracil
BRK	Propylthiouracil
TCA	Nortriptyline

The following results are tabulated from these clinical studies:

