

One Gallon (3.78-liters)

Case Packaging: 4 Each



RenalPure[®] C-116

LIQUID ACID CONCENTRATE

NOT FOR PARENTERAL USE

WARNING: This acid concentrate product is formulated to be used in conjunction with RenalPure[®] Powder Bicarbonate (CB-Series) or Sterilyte[®] Liquid Bicarbonate (CL-Series) only in a 45X three-stream artificial kidney (hemodialysis) machine.

FINAL DIALYSATE ELECTROLYTE YIELD	mg/L
Sodium	137
Calcium	2.59
Potassium	2.09
Magnesium	0.78
Acetic Acid	430
Dextrose	100
Bicarbonate	37

CONCENTRATE COMPOSITION	(Gm/L)	Density
Sodium Chloride (NaCl)	283	1.1425
Calcium Chloride (CaCl ₂)	0.24	5.84
Potassium Chloride (KCl)	0.71	0.14
Magnesium Chloride (MgCl ₂)	1.81	0.15
Acetic Acid (CH ₃ COOH)	10.8	0.04
Dextrose (C ₆ H ₁₂ O ₆)	45.0	0.24
		1.00

WARNING: For every 45 volume parts of dialysate, mix 1 volume part of this acid concentrate and 1.72 volume parts of CB-Series or CL-Series Bicarbonate Concentrate with 42.28 volume parts of water that meets AAMI quality standard for water for hemodialysis. Mix thoroughly prior to use. Test final dialysate for conductivity and pH prior to use. Refer to manufacturer for nominal conductivity of final dialysate. Avoid excessive temperature. Do not use if package is open or damaged. Federal (USA) law restricts sale of this device by or on order of a physician. Solution is non-pyrogenic.



Wiscom MI, 48393
800-448-3353
Rev. 2/18/70

17310 - 0322111
Exp. 09-22-13 CASE 349

LOT NO.

2K 2.5Ca

Manufactured in the USA by: