

DOCUMENT NUMBER	Page
DOC-053608	1 of 2
RCE DataLink™ Service Installation Checklist	REV 01

Note: Use this document as a guide to perform the service installation for DataLink with the RetCam Envision™ system.

SECTION I: CUSTOMER & RETCAM DEVICE DETAILS [COMPLETE BEFORE INSTALLATION]

Customer Information		Device Information	
Datalink Customer Information ☐ New ☐ Existing		RetCam Envision System	
☐Account # or ☐Case#		Serial Number	
Account Name		SW Version	
Address 1		RetCam Install Date	
Address 2		Total Active Devices	
City, State		DICOM License Installed?	□Yes □No
Zip Code		Device Primary Location (NICU, PICU, OR etc.)	
Country		Notes:	
Customer Time Zone			
	Customer Point of	Contacts	
Contact	Name	Email	Phone
Main Contact			
Technical/Server			
Contact			
Interface Contact			

DataLink Package			
Description Part Number			
RetCam DataLink Setup & HL7 Inbound	SVC99DATALINKINSTALL		
RetCam DataLink Subscription - 1 yr Subscription	RCDL-ANNUAL		
RetCam DataLink HL7 Outbound Interface	SVC99DATALINKHL7OUT		

Current EMR/HIS System			
EMR/HIS System HL7 Protocol			



DOCUMENT NUMBER DOC-053608	Page 2 of 2
RCE DataLink™ Service Installation Checklist	REV 01

Remote Access for Natus Support					
Remote Access Method URL User					

Project Milestones		
	Date	Notes
Server Provisioned by Customer		
Natus Access Provided		
Connectivity Testing Begins		
Integrated Testing Begins		
Go-Live		

Advanced Datalink Interface

NOTE: The following options can/may be selected but are considered ADVANCED interface options and must be forwarded to Natus Professional Services at RetCamService@Natus.com for applicable quoting of work required.

	HL7 Outbound	Final result is sent to EMR in an ORU result message
--	--------------	--

Authorized Users

First Name	Last Name	Email (Login)	Role

Roles:
ROP Coordinator - viewing, reporting, and searching
Ophthalmologist - viewing, reporting, and searching
Imager - viewing, reporting, and searching
Administrator - Add authorized users to organization, viewing, reporting, and searching

Checklist Completion (Attach completed form to Case #)					
Contact Name Signature Date					
Natus Representative					