


EnVe Ventilator

PTV Docking Station

The PTV Docking Station supports the EnVe ventilator on either the roll stand or wall mount. Docking Stations can be strategically placed throughout a hospital for power and communication. Redocking at your destination automatically starts real-time data transfer to networked medical records and provides power to the EnVe.


The Docking Station provides information in two ways: First, the VOXP communication protocol links the EnVe ventilator to monitoring and HIS systems. Second, through the PTV Smart Card, all patient trend data can be downloaded to a computer for data analysis.




The EnVe can be easily wall mounted or mobile on the roll stand.

PTM

The optional PTM provides a comprehensive graphic package that is configurable in both appearance and layout.



One touch release for easy undocking



Summary Screen Shown

- Displays 24 Hour trended patient data
- Incorporates waveforms, loops and monitors on one screen
- Data Screens display patient information and are configurable
- Flow/Volume, Volume/Pressure Loops
- Pressure, Flow and Volume Waveforms