

Product Name: 2Cool™ Zimmer Chiller Adapter

Part Number: 1500-289-00

Description: The 2Cool Zimmer Chiller Adapter provides the ability for the user to divert a portion of the cool air from the chiller directly to the patient via a second hose that is attached via the “Y” shaped adapter. As pictured in Step 1.

Directions for Use:

Step 1: Attach the locking end of the hose into the 2Cool adapter.

Step 2: Attach the single end of the adapter onto the chiller as shown below in step 2. The adapter diverts the air to the handpiece and to the patient therefore you may need to adjust your Zimmer settings based on the needs of your patient.

Step 1



Step 2

