

All in one system
for Dental

Just plug in
for immediate use!

cello



FOCUS LIGHT TYPE

SPECIFICATION	LED LIGHT TYPE	LED FOCUS LIGHT TYPE
MODEL NO.	HL-65100	HL-65100A
DIMENSION	18.89 * 17.71 * 65.35 (in)	18.89 * 17.71 * 62.99 (in)
WEIGHT	143.30(lb)	134.48(lb)
POWER SOURCE	AC 100V / 50~60Hz	
POWER CONSUMPTION	590W	
WATER/SUCTION BOTTLE	1500ml	1500ml
STORE CONDITION	32°F~104°F(TEMP.), 10%~85%(HUMIDITY)	
	LED LIGHT	LED FOCUS LIGHT
COMPONENTS	HIGH HANDPIECE, LOW HAND PIECE, SUCTION, 3-WAY SYRINGE, SCALER, TRAY, FOOT SWITCH	

cello

DENTAL ALL IN ONE SYSTEM



O:(908) 968-3553 • M:(908)397-2552
sales12@dentalfocus.biz • www.dentalfocus.biz
Sopro Tech Support:(877)422-8366 ext 3
P.O. Box 454, Neshanic Station, NJ 08853



Dental all in one system



Low-noise built-in compressor type

CELLO is built-in oilless and low-noise compressor supplies clean and dry air all the time without having to install a large compressor separately.



No need for water and waste water pipes

CELLO does not require water and waste water pipes for water supply and suction, the basic functions of dental equipment. The central vacuum for suction is not required either.



Just plug in for immediate use

CELLO is usable upon connecting the power at the desired location for use. The castors at the bottom of the unit makes CELLO movable for use any time anywhere.

SUCTION

For suctioning foreign substances including saliva generated during treatment.

HIGH SPEED HANDPIECE

Handpiece rotating at the highspeed, (400,000 RPM)

LOW SPEED HANDPIECE

Handpiece rotating at the lowspeed, (20,000 RPM, for oral prophylaxis)

SUCTION BOTTLE

For gathering foreign substances after suctioning foreign substances including saliva generated during treatment, (max 1.5L)

WHEELS

Wheels are movable freely in the treatment room and have locking device that make them easy to use.

LED LIGHT

An essential light for treatment.

LED LIGHT FOCUS TYPE ▶



POWER SWITCH & PRESSURE GAUGE

The master power switch for the product (except light for separate on/off) and tank air pressure gauge.

3-WAY SYRINGE

The syringe available to use for 3-different way syringe; Water, air and spray.

SCALER

Available to adjust water quantity for spray and scaler power when using ultrasonic scaler.

WATER BOTTLE

For the supply necessary water during treatment, (max 1.5L)

