



# EON. Smarter Body Contouring

Discover EON →

Watch Video ▶



## Providers

Discover how EON can elevate your aesthetic practice.

+ Details



## Patients

Find out if an EON treatment is right for you.

[+ Details](#)





## EON Benefits

[Learn more](#) 

---

### Elevate Your Aesthetics Practice

EON was developed by a world class R&D team with guidance from industry leading Scientific Advisory Committee members. The result is an FDA-approved

device for abdominal fat reduction evidenced through superior results.

---

## Patient Satisfaction

In clinical trials, 96% of the patients would recommend EON. There were zero non-responders, pain-free, minimal patient downtime, very low risk of side effects, and an average fat reduction of 21.6% after a single treatment (25.3% in the lower abdomen).

---



See what our **providers & patients** have to say



## EON News + Events

[View all](#) 



### EON NEWS

What is EON?

[Read more](#)



### FEATURED PRESS

EON Debuts with Body Details

[Read more](#)



EON NEWS

Unique Benefits of EON for Doctors

[Read more](#)







Enter your email to join our newsletter →

[Patients](#)

[Providers](#)

[Before & After](#)

[EON Gallery](#)

[News & Events](#)

[About](#)

[Blog](#)

[FAQ](#)

[Press](#)

[Contact Us](#)

[Events](#)

[Request Demo](#)

---

800-979-6801

[info@dominionaesthetic.com](mailto:info@dominionaesthetic.com)



---

© 2020 Dominion Aesthetic Technologies, Inc.

[Privacy Policy](#) / [Terms + Conditions](#)

Website by DIGITECH Web Design





Providers

Overview

How EON Works

Technical Features

Clinical Data

Before + After

Events



Get More Info

# Meet the Next Evolution in Aesthetics





EON is the technological benchmark in touchless robotic laser body contouring. Industry leading laser and cooling technology is coupled with years of scientific research and clinical trials. FDA cleared for abdominal fat reduction, EON provides an unsurpassed patient experience that delivers efficacy, efficiency, and ease of use that you can trust.

[Request a Demo](#) →

[Watch Video](#) ▶



# Elevate Your Aesthetic Practice

The intelligent robotic arm scans the patient while smart sensor technology maps their topography. The result is a custom treatment delivering touchless laser and cooling with the industry-best temperature differential for maximum efficacy.


While 1064nm laser energy induces apoptosis, the patient's comfort is maintained via a revolutionary cooling system. The entire abdomen is treated in a 60-minute session. The body then naturally eliminates the fat cells through the lymphatic system with best results seen at 12 weeks.





## Elegant Technology

[Get More Info](#) 



### Robotic

We've removed the inconsistency and time-consuming operator application for body contouring. Intelligent robotics provide consistency treatment after treatment with precise laser energy heat distribution. Further, robotic efficiency treats the full abdomen in 60 minutes, resulting in quicker treatments, patient value, and better provider economics.



[View EON Gallery](#) 

## CLINICAL DATA

[Download Clinical Data Summary](#) 

## Founded in Science

Our Scientific Advisory Committee methodically led study after study to ultimately bring EON's technology to the market. The team is comprised of renowned physicians and experts in medical laser aesthetics.

### Simulated Tissue Study

Determined ideal laser wavelength for fat reduction.

### Cooling Validation Study

Optimized power delivery using EON's unique scanning and cooling mechanism of delivery on

live porcine.

### Thermal Distribution Study

Validated heat distribution over the course of a treatment in a simulated phantom.

### EON Safety & Efficacy Validation Study

Multi-site clinical study confirmed EON's safety and efficacy on patients receiving a full abdominal treatment. Resulted in FDA clearance for abdominal fat reduction.

### Thermal Penetration Study


Enhanced power protocols for maximum effectiveness using thermocouples at varying depths in patients undergoing abdominoplasties.

## Before + After Treatment Results

Untouched Patient Results at 12 Weeks After a Single Treatment

BEFORE

AFTER









BEFORE

AFTER




BEFORE

AFTER



BEFORE

AFTER





[View All Before + Afters](#) →





Enter your email to join our newsletter →

[Patients](#)

[Providers](#)

[Before & After](#)

[EON Gallery](#)

[News & Events](#)

[About](#)

[Blog](#)

[FAQ](#)

[Press](#)

[Contact Us](#)

[Events](#)

[Request Demo](#)

---

800-979-6801

info@dominionaesthetic.com



© 2020 Dominion Aesthetic Technologies, Inc.

[Privacy Policy](#) / [Terms + Conditions](#)

Website by [DIGITECH Web Design](#)

---