

# products

maruho medical offers a comprehensive portfolio of orthopedic surgery products that complement shoulder, knee, and hip arthroscopies

Home » Products » Knee » SoftLink™ Button



## SoftLink™ Titanium Button

The SoftLink Button is intended for use in the fixation of bone and soft tissue in orthopedic procedures requiring ligament or tendon reconstruction. The Titanium button comes in a fixed or adjustable loop configuration. Reusable bone tunnel reamers facilitate

[contact us](#)

### ordering information

Soft Link Button

PART#	DIAMETER (MM)	LENGTH (MM)	MATERIAL	DESCRIPTION
15-11015	4	12	Titanium	non-adjustable, 15mm UHMWPE loop
15-11020	4	12	Titanium	non-adjustable, 20mm UHMWPE loop
15-11025	4	15	Titanium	non-adjustable, 25mm UHMWPE loop