

6/10/22, 2:40 PM

CHARACTERISTICS

Configuration ring

Applications medical, for CPAP therapy

Patient condition

for sleep apnea, COPD, for pneumonia, for asthma, for snoring

Options wireless

11101033

Pulse rate Max.: 250 unit

Min.: 30 unit

Oxygen level Max.: 99 %

Min.: 70 %

DESCRIPTION

Overnight Tracking Track The Trends of Your Overnight Oxygen Levels & Heart Rate.

Built-in Memory Work Independently or Together with app

Seamless Share on Sleep Analysis Share PDF/CSV Analysis & Trend Chart

Transmissive Oximetry More Sensitive on Your Body Index and More Accurate Across All Finger Size.

Long Battery Life Rechargeable Battery, Lasts 12-16 Hours

Comfortable Wearing

Weighs Only 15 grams, Competely Free Your Finger


Track Your Oxygen Level Every Second Wearing the ring, it will automatically start to track your overnight blood oxygen saturation, heart rate and body movement continuously.

Transmissive Pulse Oximetry Using intelligent SOC chip, combined with anti-interference filtering technology, O2ring provides accurate SpO2 level and Pulse Rate measurement in real time.

Real-Time Data Management


Built in memory, works stand-alone, your information is always a quick glance away with a larger display. O2ring stores 4 groups of 10-hour data. When the device is connected to the APP via Bluetooth, the data will be uploaded automatically.

VIDEO



OTHER VIATOM BY LEPU PRODUCTS

PULSE OXIMETER



See all Viatom by Lepu products

Intensive care ventilator CPAP ventilator Hand-held pulse oximeter Viatom by Lepu medical pulse oximeter Viatom by Lepu wireless pulse oximeter Home ventilator Ventilator with heated humidifier	
Viatom by Lepu wireless pulse oximeter Home ventilator Ventilator with heated humidifier	
Viatom by Lepu home care pulse oximeter Benchtop ventilator Bi-level positive-pressure ventilator OLED pulse	oximeter
Auto CPAP ventilator Battery-powered pulse oximeter Viatom by Lepu pediatric pulse oximeter	
Viatom by Lepu wrist pulse oximeter COPD ventilator	