

0 Reviews | [Write A Review](#)

Siseca Stethoscope, Classic Dual Head - for Doctors,
Nurses, Med Students, Professional Pediatric

\$25.99

Product Code: S054

Availability: In Stock

First Available: 2022-10-20

Brand: Siseca

Qty

1

[ADD TO CART](#)

PRODUCT DESCRIPTION

Siseca Stethoscope, Classic Dual Head - for Doctors, Nurses, Med Students, Professional Pediatric, Medical, Cardiology, Home Use

Diagnose and make decisions with confidence

Hear subtle changes in patient status

Ergonomic and high profile construction

Soft-sealing eartip provide optimal comfort and excellent sound occlusion

7-year warranty

Product Description

How to change the stethoscope frequency

The membrane side (high frequencies) and the funnel side (low frequencies) can be alternated by turning the stethoscope head. The stethoscope head needs to be firmly placed on the skin while listening. When measuring arterial blood pressure, the stethoscope head can be placed under the blood-pressure cuff with the membrane facing down. You can check which side is active by locating the open/closed funnel opening, or by touching the membrane with your finger while the headphones are in your ears.

PRODUCT INFORMATION

Is Discontinued By Manufacturer : No

Manufacturer : Siseca

Reviews (0)

REVIEW

Write a review

* Your Name

* Your Review

Note: HTML is not translated!

* Rating **Bad** **Good**

Continue

OUR PRODUCTS

RELATED PRODUCTS

INFORMATION

[About Us](#)

[Delivery Information](#)

[Privacy Policy](#)

[Terms & Conditions](#)

CUSTOMER SERVICE

[Returns](#)

[Site Map](#)

EXTRAS

[Brands](#)

[Gift Certificates](#)

[Affiliate](#)

[Specials](#)

[My Account](#)

[Returns](#)

MY ACCOUNT

[My Account](#)

[Order History](#)

[Wish List](#)

[Newsletter](#)

