



# InstaStain® Ethidium Bromide 100 gels 2002

**InstaStain® Ethidium Bromide, which contains only a few micrograms of ethidium bromide, is rapid and sensitive. In 2 minutes, an agarose gel is ready for visualization. Disposal of InstaStain® Ethidium Bromide is minimal compared to the volume of liquid waste generated from the standard ethidium bromide staining procedure.**

7 x 7 cm



Secure Online Ordering:

Code	Price	Qty
2002 InstaStain® Ethidium Bromide 100 gels	<b>\$89.00</b>	

[You may also place orders by phone, fax, email or regular mail](#)

(t) 1-800-338-6835 or 202-370-1500 | (f) 1-202-370-1501 | info@edvotek.com | © 2011 Edvotek Inc. All rights reserved.

### EXPERIMENTS

- Introduction to DNA
- Discovering DNA Electrophoresis
- DNA Electrophoresis (Ready-to-Load)
- Advanced DNA Applications
- Polymerase Chain Reaction (PCR)
- Introduction to RNA
- Forensic Science
- Genetic Engineering & Transformation

- Immunology
- Medical Diagnostics
- Plants
- Environmental Monitoring & Protection
- Proteins, Enzymes & Chromatography
- Advanced Placement Biology
- Reagents, Biologicals & Supplies
- Curriculum Modules

### EQUIPMENT

- DNA Electrophoresis
- Protein Electrophoresis
- Power Supplies
- Pipets & Liquid Handling
- PCR Equipment
- Microcentrifuges
- Gel Visualization
- Waterbaths
- Incubation Oven & Shakers

- Spectrophotometers
- Balances
- Meters & Timers
- LabStation™ Packages

### ABOUT OUR COMPANY

- About Us
- Contact Us
- Catalog Request