


WorkHorse® II

PERCUTANEOUS TRANSLUMINAL
ANGIOPLASTY BALLOON



**BALLOON
DIAMETER** **10** MM

**BALLOON
LENGTH** **4** CM

**CATHETER LENGTH
(Usable)** **50** CM

**RATED
BURST PRESSURE** **14** ATM

**RECOMMENDED
INTRODUCER** **7** F

**RECOMMENDED
GUIDE WIRE** **.035** IN.

**SHAFT
DIAMETER** **5** F

Catalog No. **REF** **16500444**

LOT **SAMPLE**

Use By  **2013-03**

QTY **1**

Quantite/Cantidad/Menge/
Quantita/Quantidade/Antal/Antal/
Hoeveelheid/Ποσότητα/数量

• DO NOT REUSE • STERILIZED USING ETHYLENE OXIDE • DO NOT RESTERILIZE • CAUTION
• KEEP DRY • UPPER TEMPERATURE LIMIT • KEEP AWAY FROM SUNLIGHT • DO NOT USE IF
PACKAGE IS DAMAGED • STERILE UNLESS PACKAGE IS OPENED OR DAMAGED



CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.