

Veneers

DenMat Lab was the first dental laboratory to design and fabricate thin veneers. Backed by over 30 years of experience, world-class education, and the latest digital lab technology, Lumineers produce naturally beautiful smiles with the greatest preservation of tooth structure possible. Send us a case to find out why Lumineers is the #1 patient-requested veneer.

Now Available in:

Multiple Materials

- Cerinate® II Lucite-reinforced porcelain
- IPS e.max® lithium disilicate

Multiple Preparation Techniques

- No prep
- Enamel contouring
- Traditional

Multiple Fabrication Techniques

- Pressed
- Cut-back
- Stacked

LUMINEERS

The *1 patient-requested veneer

Crowns, Bridges & Inlays/Onlays

DenMat Lab has been working with zirconia for many years. Recently, we introduced LumiZir full-contour restorations. Full-contour restorations are made from one monolithic solid piece of zirconia. We heard the request for strength without compromising esthetics. LumiZir full contour crowns are made of the strongest and most translucent zirconia on the market today.

- Strength of over 1150 Mpa flexural strength*
- High translucency
- Up to 14 units
- LumiZir crowns and bridges can be made using a LumiZir core with a porcelain or IPS e.max ZirPress®

Now Available in:

Crowns

- LumiZir
- IPS e.max
- Porcelain-fused-to-metal (PFM)
- Captek®
- Full Cast

Bridges

- LumiZir
- IPS e.max (3-unit max)
- PFM
- Maryland

Inlays/Onlays

• IPS e.max



Snap-On Smile

Many dentists regard temporization as one of the most difficult aspects of full-mouth restoration. Snap-On Smile is not only an excellent provisional, it also provides patients exceptional esthetics. It can last for years, yet is affordable enough to be temporary. For patients who are hesitant to commit to more involved restoration treatments, Snap-On Smile gives them the experience of living with better function and superior esthetics.

- Full or partial arch
- Cosmetic enhancement
- Implant provisional
- Removable partial denture
- Non-invasive

97% Excellent or Good by Clinicians Report[†]

81% of 31 CR Evaluators stated they would incorporate Snap-On Smile into their practice.

Snap-On Smile was primarily used by CR Evaluators for:

- Provisional restoration
- Short-term esthetic improvement
- Long-term esthetic improvement
- Diagnostic set-up



Data on file.

†An independent, non-profit dental education and product testing foundation: Clinicians Report,® February, 2012. For the full report contact DenMat.