Search

home

products

resources

about RMS

investors

contact us



product details > features & benefits > related products & kits > ordering information >

Site Map :: About RMS Medical Products :: Legal Disclaimer :: Privacy Policy :: Contact RMS Medical Products © 2014 RMS Medical Products. All Rights Reserved.

RES-Q-VAC Hand Held Medical Suction

performance >

- Powerful wall suction performance >600mmHg
- Works in any position will never overflow
- Essential for trach, ventilator and respiratory patients
- Extremely durable and reliable ready when you need it

cost & value >

- RES-Q-VAC® is a fraction of the cost of electric counterparts
- No maintenance costs involved with checking suction levels, unlike electric units

safety >

- The available 0.22 micron (μ) Full Stop Protection® (FSP) Filter protects against infectious diseases including Hepatitis, Tuberculosis, HIV and H1N1
 FSP protects the pump for safe reuse
- · Sterile and non-sterile configurations
- RES-Q-VAC® has full tactile feedback and whistle-tip catheters for patient safety

convenience >

- Portable hand-powered suction
- · No batteries or electric required
- Ideal for crash carts, ER, patient transport, patients in isolation, cafeterias
- · Full line of push connector, whistle tipped catheters

<u>Literature Instructions Photos Videos</u>

fin t




					Search
home	products	resources	about RMS	investors	contact us



product details >
 features & benefits >
related products & kits >
 ordering information >

RES-Q-VAC Hand Held Medical Suction



To find out more, please call us toll-free at 800-624-9600

Part#	Description	Part#	Description
Typical Ho	ospital Configuration	S	
P2F1B1AFS	RES-Q-VAC™ Pump	R1F1B1AFS (replacement)	FSP Canister
	2 FSP Canisters		Ridgid Yankauer
	Ridgid Yankauer		Soft Yankauer
	Soft Yankauer		3.5"/14" 8fr Catheters
	3.5"/14" 8fr Catheters		3.5"/14" 10fr Catheters
4 w 100 c 1 m 100 c 100	3.5"/14" 10fr Catheters		3.5"/14" 14fr Catheters
	3.5"/14" 14fr Catheters		Sterile
	Carry Bag		
	Sterile		
P1F1B1AFS	RES-Q-VAC™ Pump	R1F1B1AFS (replacement)	FSP Canister
And the second s	FSP Canister		Ridgid Yankauer
	Ridgid Yankauer		Soft Yankauer
	Soft Yankauer		3.5"/14" 8fr Catheters
	3.5"/14" 8fr Catheters		3.5"/14" 10fr Catheters
TOTAL CONTRACTOR OF THE PARTY O	3.5"/14" 10fr Catheters		3.5"/14" 14fr Catheters
androgen and an arrange and an arrange and arrange and arrange and arrange and arrange and arrange and arrange	3.5"/14" 14fr Catheters		Sterile
	Carry Bag		
	Sterile		
P1F1S1C1S	RES-Q-VAC™ Pump	R1F1S1C1S (replacement)	FSP Canister
	FSP Canister		Soft Yankauer
***************************************	Soft Yankauer	3.5"/14"	3.5"/14" 10fr Catheters
	3.5"/14" 10fr Catheters		14" 14fr Catheters
in a transfer of the second	14" 14fr Catheters		Sterile
	Саггу Вад		
	Sterile		
Typical Ho	ome Care Configurat	ions (Ventil	ator Patient)
P1N001145	RES-Q-VAC™ Pump	R1N00114S (replacement)	Non-FSP Canister
	Non-FPS Canisters		3.5"/14" 14fr Catheters
	3.5"/14" 14fr Catheters		Sterile
	Carry Bag		
	Sterile		

Typical Ho	ome Care Configurati	ons (Trache	eostomy Patient)
P1N001815	RES-Q-VAC™ Pump	R1N00181S (replacement)	Non-FSP Canister
,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	Non-FPS Canisters		3.5"/14" 8fr Catheters
	3.5"/14" 8fr Catheters		3.5"/14" 10fr Catheters
	3.5"/14" 10fr Catheters		Sterile
	Саггу Вад		
***************************************	Sterile		
Typical Fi	rst Responder Config	urations	
P1F1S110N	RES-Q-VAC™ Pump	R1F1S110N (replacement)	FSP Canister
	FPS Canisters		Soft Yankauer
	Soft Yankauer		3.5"/14" 10fr Catheters
	3.5"/14" 10fr Catheters		
	Carry Bag		
P1F1R110N	RES-Q-VAC™ Pump	R1F1R110N (replacement)	FSP Canister
***************************************	FPS Canisters		Soft Yankauer
	Soft Yankauer		3.5"/14" 10fr Catheters
	3.5"/14" 10fr Catheters		
Typical N	ursing Home Configu	rations	
P1F1S114S	RES-Q-VAC™ Pump	R1F1S114S (replacement)	FSP Canister
	FPS Canisters		Soft Yankauer
	Soft Yankauer		3.5"/14" 14fr Catheters
	3.5"/14" 14fr Catheters		Sterile
	Carry Bag		
	Sterile		

<u>Literature</u> <u>Instructions</u> <u>Photos</u> <u>Videos</u>

Site Map :: About RMS Medical Products :: Legal Disclaimer :: Privacy Policy :: Contact RMS Medical Products © 2014 RMS Medical Products. All Rights Reserved.

Search

home

products

resources

about RMS

investors

contact us

RMS resources



literature >

instruction >

product photos >

videos >

FREEDOM60" & **Precision Flow** Rate Tubing™

HIgH-Flo Subcutaneous Safety Needle Sets"

RES-Q-VAC® Hand Held **Medical Suction**



RES-Q-VAC Photos >

- RES-Q-VAC Complete Package
- **RES-Q-VAC Replacement Package**
- **DENTAL-E-VAC Complete Package**
- **Snapping on Canister**
- **Attaching Catheter**
- Squeezing RES-Q-VAC Pump
- RES-Q-VAC Pump Handle
- Canister Replacement Kit
- **Full Stop Protection Cap**
- **RES-Q-VAC Canister** Patient Label

- Sampler Bottle RES-Q-VAC with LED and Rigid Yankauer
- RES-Q-VAC with LED On
- RES-Q-VAC LED
- **RES-Q-VAC First Responder Use**
- Physician Using RES-Q-VAC on Baby
- RES-Q-VAC in Use

training >







Site Map :: About RMS Medical Products :: Legal Disclaimer :: Privacy Policy :: Contact RMS Medical Products © 2014 RMS Medical Products. All Rights Reserved.

