

100g (3.5 oz) / Bag

Case Packaging: 2 Each



RenalPure®

POWDER BICARBONATE

CB-25

NOT FOR PARENTERAL USE

WARNING: This bicarbonate concentrate product is formulated to be used in conjunction with DeSorb® Dry Acid (DC-Series) or RenalPure® Liquid Acid (C-Series) only in a 45X three-chamber artificial kidney (hemodialysis) machine.

ELECTROLYTE YIELD (mixed 1:1)	
Sodium	37.0 mEq/L
Bicarbonate	37.0 mEq/L

CONCENTRATE COMPOSITION	
Per Bag (Grams)	(Gentl) Concentrate
Sodium Bicarbonate (NaHCO ₃)	7892
	Dilute
	125.16
	3.11

WARNING: This carton will produce 50 gallons (189.25 liters) of liquid bicarbonate concentrate when entire contents are properly mixed. Use water that meets AAMI quality standards for water for hemodialysis. Read mixing instructions prior to use. Once mixed, product should be used within 24 hours. Bacterial growth can occur in liquid bicarbonate concentrate. Residual liquid bicarbonate concentrate and poorly cleaned mixing tanks and systems will support bacterial growth. Test final dialysate for conductivity and pH prior to use. Refer to manufacturer for nominal conductivity of final dialysate. Avoid excessive temperature. Do not use if package is open or damaged. Federal (USA) law restricts sale of this device by or on order of a physician.



ROCKWELL

MEDICAL

Wixom MI 48393
800-449-3353
Rev. 1/1/10

17331 - 0324111

818

Manufactured
in the USA by

LOT NO.