

# Endovascular Catalog - U.S.

## SLEEK® OTW PTA Dilatation Catheter

Designed for infrapopiteal vessels. Delivers precision, speed and control.

### .014" Over-The-Wire

- Coaxial
- High rated burst pressure (up to 16atm)
- Long Balloon lengths (1.5-28cm)

Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Shaft Length (150cm)
1.25	15	14	4	4261201X
1.5	15	14	4	4261501X
1.5	20	16	4	4261502X
1.5	40	16	4	4261504X
1.5	80	16	4	4261508X
1.5	100	16	4	4261510X
1.5	120	16	4	4261512X
2	20	16	4	4262002X
2	40	16	4	4262004X
2	80	15	4	4262008X
2	100	15	4	4262010X
2	120	15	4	4262012X
2	150	15	4	4262015X
2	220	15	4	4262022X
2	280	15	4	4262028X
2.5	20	16	4	4262502X
2.5	40	16	4	4262504X
2.5	80	15	4	4262508X
2.5	100	15	4	4262510X
2.5	120	15	4	4262512X
2.5	150	15	4	4262515X
2.5	220	15	4	4262522X
2.5	280	15	4	4262528X

Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Shaft Length (150cm)
3	20	16	4	4263002X
3	40	16	4	4263004X
3	80	15	4	4263008X
3	100	15	4	4263010X
3	120	15	4	4263012X
3	150	15	4	4263015X
3	220	15	4	4263022X
3	280	15	4	4263028X
3.5	20	16	4	4263502X
3.5	40	16	4	4263504X
3.5	80	15	4	4263508X
3.5	100	15	4	4263510X
3.5	120	15	4	4263512X
3.5	150	15	4	4263515X
3.5	220	15	4	4263522X
3.5	280	15	4	4263528X
4	20	16	4	4264002X
4	40	16	4	4264004X
4	80	15	4	4264008X
4	100	15	4	4264010X
4	120	15	4	4264012X
4	150	15	4	4264015X
4	220	15	4	4264022X
4	280	15	4	4264028X
5	20	14	4	4265002X
5	40	14	4	4265004X
5	80	13	4	4265008X
5	100	13	4	4265010X
5	150	10	4	4265015X
5	220	10	4	4265022X
5	280	10	4	4265028X

## SLEEK® RX PTA Dilatation Catheter

Designed for infrapopiteal vessels. Delivers precision, speed and control.

### .014" Rapid Exchange Balloon

- Dual-lumen rapid exchange shaft
- High rated burst pressure (16-15atm)
- Long balloon lengths (4-22cm)

Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Shaft Length (150cm)	Shaft Length (155cm)
2	40	16	4	4252004X	
2	80	15	4	4252008X	
2	100	15	4	4252010X	
2	120	15	4	4252012X	
2	150	15	4		4252015X
2	220	15	4		4252022X
2.5	40	16	4	4252504X	
2.5	80	15	4	4252508X	
2.5	100	15	4	4252510X	
2.5	120	15	4	4252512X	
2.5	150	15	4		4252515X
2.5	220	15	4		4252522X
3	40	16	4	4253004X	
3	80	15	4	4253008X	
3	100	15	4	4253010X	

Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Shaft Length (150cm)	Shaft Length (155cm)
3	120	15	4	4253012X	
3	150	15	4	4253015X	
3	220	15	4	4253022X	
3.5	40	16	4	4253504X	
3.5	80	15	4	4253508X	
3.5	100	15	4	4253510X	
3.5	120	15	4	4253512X	
3.5	150	15	4		4253515X
3.5	220	15	4		4253522X
4	40	16	4	4254004X	
4	80	15	4	4254008X	
4	100	15	4	4254010X	
4	120	15	4	4254012X	
4	150	15	4		4254015X
4	220	15	4		4254022X

For information on indications, contraindications, warnings, and precautions, see page 114.

## U.S. Headquarters

1820 McCarthy Blvd.  
Milpitas CA, 95035  
Phone: 408.273.3700

## Customer Service Center

Phone: 1.800.327.7714 (toll-free, continental U.S.)  
Fax: 800.997.1122

Attn: Customer Contact Center  
14201 Northwest 60th Ave  
Miami Lakes, FL 33014

## Technical Information

1.800.327.7714 (toll-free, continental US)

All orders should be directed to the Customer Service Center by fax, EDI, or telephone. The Customer Service Center is open weekdays (excluding US national holidays) between the hours of 7:30 AM and 9:00 PM EST/EDT. Emergency order service is also available 24 hours a day.

## Terms

Net 30 days. Normal delivery second day A.R.O., F.O.B. shipping point with freight prepaid. Standing order shipments are F.O.B. shipping point with freight prepaid. Prices are subject to change without notice.

## Disclaimer of warranty and limitation of remedy

There is no express or implied warranty, including any implied warranty of merchantability of fitness for a particular purpose, including the products described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

## Standing Purchase Orders

Cordis encourages standing purchase orders. Standing purchase orders are agreements to purchase specific quantities of product for delivery on a specific, predetermined schedule. To establish a standing purchase order, telephone the Customer Contact Center at 1.800.327.7714 option 1, or contact your local Cordis Sales Representative.

## Intended Product Usage/Storage

- **Usage:** Most products in this catalog are supplied sterile and nonpyrogenic. Do not use any products if their sterile package is damaged. Discard catheters and single-use accessories after one procedure. All parts are extremely difficult to clean. DO NOT REUSE OR RESHAPE. DO NOT AUTOCLAVE SINGLE-USE PRODUCTS. Structural integrity and/or function may be impaired through reuse, reshaping, and cleaning. ACCORDINGLY, CORDIS CORPORATION WILL NOT BE RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM REUSE

OF CATHETERS AND ACCESSORIES. Prior to use, refer to the instructions accompanying the product.

- **Storage:** Store products in a cool, dark, dry place. Use sterile products prior to the "Use By" date. Do not expose to organic solvents.
- **Caution:** Federal (USA) law restricts these devices to sale by or on the order of a physician.

## Return Policy

Products on consignment are not returnable for credit under this policy. To return any product on consignment, you must contact your local Sales Representative for instructions.

## Return Authorization

To return other Cordis products you must contact the Customer Service Center at 1.800.327.7714 option 2 or your local sales representative for a return goods authorization. You must also provide information about the acquisition method (e.g. purchase, consignment, or evaluation) for the product being returned.

## Credits

Full Credit of the invoice price will be issued on products returned within 180 days from the date of invoice. Full Credit of the invoice price, less a 10% restocking charge, will be issued on products returned later than 180 days from the date of invoice with the exception of EXOSEAL® Vascular Closure Device.

Partial credit will be issued for product returned in quantities less than the full 10 pack for EXOSEAL® Vascular Closure Device when returned within 365 days.

For customer service, call 1.800.327.7714. To learn more, visit [cordis.com](http://cordis.com).  
CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.  
See package insert for full product information. All claim data is on file at Cordis.  
© 2018 Cardinal Health. All Rights Reserved.  
100510144 3/2018