



Library

Products

Drug Search

Training

Contact Us

General

Drug Testing

Alcohol Screen

Validity Screens

Clinical Controls

Specials

**iCup<sup>A.D.</sup>**  
Drug Screen<sup>®</sup>

[Instructions](#) [Product Brochure](#) [Information Library](#)



Picture 1 | Picture 2 | Picture 3 | 360°

**Is it possible to cheat the iCup<sup>®</sup> drug screen by substituting or adulterating the specimen?**

No. The iCup<sup>™</sup> A.D. contains a specimen validity test or SVT strip that verifies the validity of a specimen. The SVT strip checks for oxidants, specific gravity and pH in the specimen. For example, if someone tries to substitute water for urine, the specific gravity indicator will show an abnormal reading.

[more frequently asked questions](#) ▶

Detects the presence of drugs.  
It's drug testing efficiency you can count on.

**Integrated testing solution:**

- **iCup<sup>®</sup> with temperature strip**
- **Two-part result form available**
- **Security Seal**
- **Photocopy Template**

**Benefits:**

- **Room temperature storage**
- **Built-in validity test**
- **Detects 3-10 drugs**
- **24-month shelf life**
- **Simple procedure**
- **Fast results**

Drug Name:	Abbr:	Cutoff:
Amphetamine	AMP	1000ng/mL
Barbiturates	BAR	300ng/mL
Benzodiazepines	BZO	300ng/mL
Cocaine	COC	300ng/mL
Ecstasy	MDMA	500ng/mL
Marijuana	THC	50ng/mL
Methadone	MTD	300ng/mL
Methamphetamine	mAMP	1000ng/mL
Morphine	MOR	300ng/mL

Opiates	OPI	2000ng/mL
Oxycodone	OXY	100ng/mL
Phencyclidine	PCP	25ng/mL
Propoxyphene	PPX	300ng/mL
Tricyclic Antidepressants	TCA	1000ng/mL

Office Hours: 8:30am to 6:00pm EST

---

[Home](#) | [Library](#) | [Products](#) | [Drug Search](#) | [Training](#) | [Contact Us](#) | [General](#) | [Site Map](#)  
© 2008 Instant Technologies, Inc. All Rights Reserved