







CoolTouch II™, the premier system for CoolTouch® Wrinkle Treatments, is designed to offer effective skin rejuvenation that is safe for all skin types. When used with the TQ10 cooling accessory, the CoolTouch II sequences pulsed 1320nm light with powerful pulsed cooling. Absorbed by water, the pulsed 1320 light heats a uniform layer of upper dermis to therapeutic levels, while the pulsed cooling effectively protects the epidermis. A thermal sensor measures the immediate response and assists the operator in delivering both safe and effective CoolTouch wrinkle treatments.

Benefits of the CoolTouch II System:

- No downtime for patients
- All skin types can be treated



CoolTouch II Specifications

Wavelength

1320nm

Fluence

Up to 24 J/cm²

Beam Delivery S-10 CoolTouch Fiber Delivery Device

Dimensions

31" high, 18" wide, 21" deep

(79cm x 46cm x 53cm)

Weight

160 lbs (73KG)

Utilities

115 VAC, 50/60 Hz, 20A, Single Phase

230 VAC, 50/60 Hz, 15A, Single Phase

Accessories

TQ10

Cooling handpiece with pre-cooling

and Thermal Quenching™ with 10mm spot size

DCD

Cooling handpiece







*CoolTouch is the 1st infrared laser system cleared for peri-orbital wrinkle treatment.

© 2001 CoolTouch Corporation, 9085 Foothills Blvd., Roseville, CA 95747 ph: 916.677.1975 • fx: 916.677.1976 • toll-free: 877.858.2665 • www.cooltouch.com P/N 7110-0002 RvA



From the Company
that Introduced
Noninvasive
Skin Rejuvenation

