



CEREC Applications in TMD and Sleep Apnea Treatment Planning

Presented by

Dr. Jason McMillan, DMD

mint

Mint Dental Works

Portland, Oregon



SICAT FUNCTION Advantages

- Fully digital treatment for TMD
- Minimal patient appointments
- Minimal “doctor time”
- Precise recordings of all movements of the jaw
- Real condyle-fossa relationship

mint

SICAT FUNCTION Advantages

- Therapeutic appliance requiring minimal adjustments

SICAT OPTIMOTION
Appliance



mint

SICAT FUNCTION Workflow

- Optical surface scan of upper and lower arch (full arch scans in CEREC)




- No messy impressions
- No models to pour
- No model inaccuracies

mint

SICAT FUNCTION Workflow


- SICAT Fusion Bite plate is lined with bite registration material for upper and lower jaw. The SICAT Fusion Bite is worn during the CBCT scan and SICAT JMT recording of the Fusion Bite



mint

SICAT FUNCTION Workflow

- CBCT scan taken with patient wearing SICAT Fusion Bite



mint