
Vouchers


TrueStim TREAT · RECOVER · PREVENT LogIn



PURCHASE Why TrueStim? News & Articles Testimonials

Customer Support Careers Distributor Opportunities






TENS

Pain signals travel through nerve receptors to the brain through the nervous system. We deliver bio-electric pulses and rhythms that mimic the body's natural signals so you feel our technology instead of the pain sensations. This has a lasting affect even after using the TrueStim.

EMS

Electric muscle stimulation trains your muscles in a way that traditional workouts alone can't. EMS, also known as neuromuscular electric stimulation, muscle stim, or e-stim, is delivered through electronic pulses to your nerve fibers in order to create involuntary muscle contractions. Unlike TENS, EMS waveforms activate fast twitch and slow twitch muscle fibers allowing the muscles to be stimulated to the point of exhaustion (relaxation) thus healing them to become stronger. This is similar to weight training. Muscles support joints and ligaments creating better support.





Microcurrent

Microcurrent treatments concentrate on pain relief and speeding healing and recovery. It is commonly used by professional and performance athletes with acute pain and/or muscle tenderness as it is drug-free and non-invasive, thus avoiding testing and recovery issues. It is also used as a cosmetic treatment. Mode 11 on the wireless touch is the microcurrent setting. Best used with the ear clip attachment.

Relieve Your Pain Contact Increase Muscle Performance Careers Instructions **Distributor Opportunities Jacksonville Jaguars Register Device Warranty Return Policy**



© TrueStim, Inc. 2017