

Equipment Refurbishing




Celebrating 30 Years of Refurbishing Excellence

Since 1984, Avante Medical Surgical has prided itself on the superior quality of its refurbishing process. We have decades of experience refurbishing monitors, anesthesia machines, defibrillators, electrosurgical equipment, surgical lights, and much more.

Our Philosophy


When it comes to finding a refurbishing partner, it is important to choose a company that can keep up with the demanding nature of the business. Avante Medical Surgical takes the extra steps to be the supplier you can trust.



Which Refurb Level is Right for Your Facility?

Reprocessed Equipment devices are checked to see if they meet original equipment manufacturer standards. After that, they're cleaned, disinfected, sterilized and shipped to your door.

Cosmetically Refurbished Equipment devices are tested to meet OEM specifications. but biomedical engineers paint and complete other surface



Which Refurb Level is Right for Your Facility?

Reprocessed Equipment devices are checked to see if they meet original equipment manufacturer standards. After that, they're cleaned, disinfected, sterilized and shipped to your door.

Cosmetically Refurbished Equipment devices are tested to meet OEM specifications, but biomedical engineers paint and complete other surface modifications before it leaves the warehouse.

Parts Refurbished Equipment the device didn't meet OEM specifications, biomedics will replace the necessary parts then cosmetically refurbish the equipment.

Fully Refurbished Equipment the ultimate overhaul of a device, biomedics replace major parts and give the equipment a full facelift to meet or exceed OEM specifications.

ANESTHESIA MACHINES

SURGICAL TABLES

Respiratory Ventilators Our Refurbishing Process

- 1 Grounding resistance is checked, current leakage is tested.
- 2 Housing is inspected and cleaned or replaced.
- 3 [Blank]
- 4 [Blank]
- 5 [Blank]
- 6 [Blank]
- 7 [Blank]
- 8 All alarms and audible signals are tested and verified.
- 9 Performance and maintenance tests are performed.

SURGICAL LIGHTS

2 Controls and switches are inspected and replaced, if necessary.

5 Performance and illumination is tested, repaired if necessary.

3 Existing bulbs and fiber optic cables are replaced.

6 The unit receives a new label and is shipped to the customer.

ESUS

RESPIRATORY VENTILATORS

DEFIBRILLATORS

Looking for something in particular?

Call: (800) 462-8195

Our Expertise

- Medical Surgical
Patient Monitoring
Ultrasound

Products & Services

- Rental Services
Avante International
Global Outreach & Development

Account

Sign up / Log in

Help

- Customer Service
Technical Support
Video Training

Avante

- About
Commitment to Quality
Careers