

This site uses cookies to provide you a better user experience. You hereby consent to our cookies if you continue to use our website. More information on our [Cookie Policy](#), [Privacy Policy](#) and [Terms of Service](#)



English



MEDICAL SPECIALTIES

PRODUCTS

EVENTS & EDUCATION

CUSTOMER SERVICE

ABOUT US

[Home](#) > [Products](#) > [Orthopedics](#) > [Hip](#)

PressFT™ and GENESYS™ PressFT™ Suture Anchors



Clinical Resources

[Videos and Documents](#)

[Instructions for Use](#)

Need More Information? Contact Us

Connect with your CONMED
Sales Representative

Product Overview

For surgeons who need strong fixation in a small footprint, the PressFT™ Suture Anchors combine small size with exceptional strength. The implant comes single or double-loaded with Hi-Fi® suture and is available in PEEK or GENESYS® biocomposite material.

PressFT can be used for Hip and Shoulder Labral repair as well as MPFL Reconstruction.

Videos & Documents



Specifications and Ordering Information



Customer Service



CONMED is committed to providing our customers with best-in-class customer service. Below are some quick links to CONMED product information that is easily accessible online. For additional assistance or to speak to a CONMED Customer Service Representative to order these CONMED products [click here](#).