



# MONTHLY CONDITIONING SOLUTION

For Fluorosilicone  
Acrylate RGP Lenses  
Coated With

**tangible**  
HYDRA-PEG®



Restores the Hydra-PEG Coating

Sterile .34 FL OZ (10 mL)

1 Treatment | 1 Lens Case

FEDERAL LAW RESTRICTS THIS  
DEVICE TO SALE BY OR ON THE ORDER  
OF A LICENSED PRACTITIONER  
This device is intended for use by a Tangible Science  
licensed practitioner only. © 2019 Tangible Science, Inc. All rights reserved.

Tangible Boost  
Subject 33

Tangible Science  
Redwood City, CA 94063  
Caution: Investigational  
United States Law for  
Investigational Use Only

*Hey Tangible Boost, thanks for everything...you're the best!  
xxxxx  
Left Eye & Right Eye  
(Your Biggest Fan!)*