

Q

PRODUCTS

EQUINOX LIBRARY AMPLIFICATION KIT

Setting a new standard in NGS Library Amplification

Kit contents

Equinox Library Amplification Mix (2X), P5/P7 Primer Mix (10X, optional)

1 of 3 8/23/2022, 10:39 AM



Q

Storage	Long term storage: -20°C ± 5°C; Short term use: 4°C for up to 2 weeks
Shelf life	12 months
Stability	25 freeze-thaw cycles



WATCHMAKER GENOMICS 5744 CENTRAL AVE. SUITE 100 BOULDER, CO 80301

720-543-2174

TECHNOLOGY

PRODUCTS

CUSTOM SOLUTIONS

COMPANY

CAREERS

NEWS AND EVENTS

2 of 3 8/23/2022, 10:39 AM



C

GENERAL REQUEST

REQUEST A QUOTE

TECHNICAL SUPPORT

CUSTOM SOLUTIONS

CERTIFICATE OF ANALYSIS

NEW CUSTOMER ACCOUNT

Join our newsletter

Enter your work email				
I'm not a robot	reCAPTCHA Privacy - Terms			

SUBMIT



Privacy Policy Terms of Use

3 of 3 8/23/2022, 10:39 AM