

search...

[Versión en Español]

Home

Advocate

Advocate Duo

Redi-Code

Other Advocate Products

Web Links

Sales/Service

To assure the influent even of eathboulds only assure the influent even of eathboulds only assure the influent even of eathboulds only assure the influent even of eathboulds our assure the influent even of eathboulds our

मान्यामातु ५०० नमान्य, कार्य मधीतु वार्यानस्थानस्य मान् वहार वा मान् नद्यावामानस्य ५०० महरूरः औरत्यः अस्योऽ भारता परा वृद्धास्त्रः जा महरू व्यवद्यव्यक्षः



Strips and Control



## **Test Strips & Control Solution**

Advocate Quick Draw strips make the test easy and accurate

- Easy to handle
- Flip-Top bottle design: Opens, closes and handels easily.
- Patented bottle design with built-in desicant: Keeps strips fresh after opening.
- Capillary Action Technology and Visual Confirmation: You can see how quick a Advocate strip automatically draws the right amount of blood for an accurate reading.



## Main Menu

- Home
- Advocate
- Advocate Duo
- Redi-Code
- Other Advocate Products
- Web Links
- Sales/Service
- Online Log Book
- BP Warranty Card

Who's Online

We have 8 guests online

Site Stats

Content View Hits: 115541



## **Control Solution**

Control Solution is used to prove the accuracy of the Advocate meter and test strips. Control solution has three functions:

- 1. To practice doing a test without using blood.
- 2. To ensure that your meter and tests strips are working together properly.
- To verify that you are doing the testing procedure properly.



Strips/Control

Software/Cable

**Lancing Products** 

Heating Pad

Policies

## Strips and Control Solution

Advocate Quick Draw strips make the test easy and accurate. Control Solution is used to prove the accuracy of the Advocate meter and test strips. Learn more...

© 2008 Diabetic Supply Of Suncoast