

*** User:soestreich ***

#	Total Marks	Dead Marks	Live Viewed Docs	Live Viewed Images	Status/ Search Duration	Search
01	1	0	1	1	0:01	88455594[SN]
02	60	0	60	54	0:01	*fl{v}{ "ckqxs"1:2}\$1pr*[bi,ti] not dead[ld]
03	0	0	0	0	0:01	*phl{v}{ "ckqxs"1:2}\$1pr*[bi,ti] not dead[ld]
04	10910	N/A	0	0	0:01	*fle{ "ckqxs"1:2}*[bi,ti] not dead[ld]
05	75383	N/A	0	0	0:02	*pro*[bi,ti] not dead[ld]
06	203	N/A	0	0	0:01	4 and 5
07	15522	N/A	0	0	0:01	*fl{v}{ "ckqxs"1:2}*[bi,ti] not dead[ld]
08	258	0	258	235	0:01	5 and 7
09	23	0	23	21	0:01	*phle{ "ckqxs"1:2}*[bi,ti] not dead[ld]
10	4117	N/A	0	0	0:01	*flex*[bi,ti] and ("009")[cc] not dead[ld]
11	2461	N/A	0	0	0:01	*flex*[bi,ti] and ("007")[cc] not dead[ld]
12	1837	N/A	0	0	0:02	*flex*[bi,ti] and ("009" a b "200")[ic] not dead[ld]
13	882	N/A	0	0	0:02	*flex*[bi,ti] and ("007" a b "200")[ic] not dead[ld]

Session started 8/29/2019 4:37:45 PM

Session finished 8/29/2019 4:55:28 PM

Total search duration 0 minutes 16 seconds

Session duration 17 minutes 43 seconds

Default NEAR limit=1ADJ limit=1

Sent to TICRS as Serial Number: 88455594