



SCHILLER

CARDIOVT MS-2010

15:48

Resting ECG

120 /min

HR 120 /min  
P 88 ms  
PR 136 ms

RRS 500 ms  
P 88 ms  
FR 136 ms

QRS 80 ms  
QT 270 ms  
QTC 386 ms

P 29°  
QRS 46°  
T 40°

R L F  
C1 C2 C3 C4 C5 C6

5 mm/2V 50 mm/s

0.15 - 35 Hz



Start

Manual

Copy

Patient

Menu

0.05 - 150 Hz F50  
SCHILLER SWITZERLAND  
Art. No. 2.18.100  
CE 0122 9631