




JET FREE FLOW™ CATHETERS



 Download brochure

[Ordering Information](#)

Jet Free Flow™'s Cross Section



> [tip illustration...](#)



Catheters

Contact

Ordering

FEATURES & BENEFITS

Soft Tip Stylet (Short Term) Integral soft tip reduces the risk of vessel damage.

Crack Resistant Luers The Ergonomic design ensures greater resistance to cracking or separation.

Implantable Cuff (Long Term) Implantable cuff allows for secure anchoring.

Priming Volume Priming volume designations clearly printed on catheter ID tags for ease of catheter maintenance.




Rotating Suture Wing Rotating suture wing provides secure external anchoring.

PRE-LOADED STYLET Allows standard or modified Seldinger introduction and once removed, renders the venous lumen totally open, unlike tapered tipped catheters.

NO SIDE HOLES Since venous and arterial lumens are totally open, there is no need for side holes to increase efficiency. Side hole fibrin formation, a major factor in catheter thrombus, is eliminated.

STRAIGHT AND PRE CURVED CONFIGURATIONS Free Flow™ Catheters available in straight and pre curved configurations for both long and short term patients.

BACK CURVED CATHETER BODY (Short Term) Permits internal jugular placement



supported by the clavicle. This curve helps reduce infection rate, helps eliminate kinking and positions the catheter terminations on the chest wall (unlike traditional “neckline” catheters).




Unique Tip Design > Designed for maximum flow rates & efficient clearances.

TRUE 360° FLOW

- > Utilizes proven coaxial, tube within tube, design
- > Reduces arterial insufficiencies
- > No side holes reduces the loss of locking fluid

OPEN TIPS Unlike standard tapered Seldinger Catheters, the arterial and venous lumens remain open for unobstructed uptake and return. Pressure/flow ratios are maximized for efficient clearances.





ORDERING INFORMATION

Jet Free Flow™



Catheters

Contact

Ordering

PLACE ORDER

On-line Order Form!

Download brochure

• [Catheter Information](#) •

SETS (EACH)	DESCRIPTION
Long Term	
JFFL1524S	15.5FR x 24cm Straight
JFFL1528S	15.5FR x 28cm Straight
JFFL1532S	15.5FR x 32cm Straight
JFFL1536S	15.5FR x 36cm Straight
JFFL1524PC	15.5FR x 24cm Pre-Curved
JFFL1528PC	15.5FR x 28cm Pre-Curved
JFFL1532PC	15.5FR x 32cm Pre-Curved
JFFL1536PC	15.5FR x 36cm Pre-Curved
Short Term	
JFFS1515S	15.5FR x 15cm Straight
JFFS1520S	15.5FR x 20cm Straight
JFFS1524S	15.5FR x 24cm Straight
JFFS1528S	15.5FR x 28cm Straight
JFFS1512IJ	15.5FR x 12cm Pre-Curved
JFFS1515IJ	15.5FR x 15cm Pre-Curved
JFFS1520IJ	15.5FR x 20cm Pre-Curved

Long Term

Sets Contain:

- (1) Catheter with Stylet
- (1) New Grasping Tunneling Tool
- (2) Vessel Dilators
- (1) Tearaway Sheath
- (1) J-Flex Guidewire
- (1) Introducer Needle
- (2) Injection Caps

Short Term

Sets Contain:

- (1) Catheter with Stylet
- (3) Vessel Dilator
- (1) J-Flex Guidewire
- (1) Introducer Needle
- (2) Injection Caps

For help with placing an order, send an email to:

customerservice@jetmedical.net

HOW TO **PRODUCTS** **CUSTOMER**

Enter Your Information Below

Name_

Company_

Phone_

***Fax_**

Email_

Address_

City_

State_

Zip_

Country_

***Shipping_**

Comments:

**not required*

SUBMIT FORM

(You will receive order confirmation).