## CLASS 21



#### **ELECTROSURGICAL**

#### **KICK BUCKET**



# WALLACH

## Disposable, Extra-rigid Loop and Ball Electrodes

Wallach offers a wide variety of sterile, single-use electrodes, all made in the U.S.A. and fashioned out of high-quality tungsten wire.

McKesson #	Description	Packaging	Mir#
317967	3mm Bali	5/box	909038
221813	5mm Ball	5/box	909003
475743	5mm Ball (6cm Shaft)	5/box	909136
449031	10mm x 7mm Loop	5/box	909133
513554	15mm x 5mm Loop	5/box	909134
294068	15mm x 8mm Loop	5/box	909011
221553	10mm x 10mm Loop	5/box	909007
222759	20mm x 8mm Loop	5/box	909013
557540	0.8mm Dia, X 16mm	5/box	909135
449030	10mm x 10mm Square	5/box	909133
520536	0.8mm Dia, X 25.4mm Large Needle	5/box	909005
209135	10mm x 8mm Square	5/box	20.00
435993	20mm x 10mm Loop	5/box	909017
221554	20mm x 15mm Loop	5/box	909132 9 <b>0</b> 9009



### **M**SKESSON

#### entrust™ Performance Stainless Steel Kick **Bucket with Frame**

360° flat bumper on frame. Heavy-duty, non-corrosive stainless steel construction. Not made with natural rubber latex.

McKesson #	Description	Packaging	Mtr #
407338	Kick Bucket	each	81-48711

#### **LAP SPONGES**



#### **INFILTRATION PUMPS**





#### Klein Infiltration Pump

An essential device designed specifically to expedite the delivery of tumescent anesthesia. Delivers high volumes of tumescent anesthesia, therefore maximizing patient comfort and reducing physician infiltration time. The only infiltration pump with a non-conductive roller head, a feature that ensures safety during procedures. The pumps are stackable for conserving space. Includes instructional DVD. 1 year warranty.

McKesson #	Description Klein Infiltration Pump	Packaging	Mfr#
887771		each	KJP II
Infiltration Tu			
McKesson #	Description Single Spike, 15.5 ft Single Spike, 15.5 ft Single Spike, No-Drip, 15.5 ft	Packaging	Mfr #
863937		10/box	ITS-10
863938		20/box	ITS-20
882686		10/box	CVT-10
957266	Single Spike, No-Drip, 15.5 ft	20/box	CVT-20
863936	Double Spike, 15.5 ft	20/box	ITD-20

#### **MCKESSON**

#### Medi-Pak™ Lap Sponges, Prewashed

Description

100% cotton, X-Ray detectable, Sterile. Single use. Not made with natural rubber latex.

Packaging

#### Soft Pack McKesson #

281500	12" x 12"	5/pack, 40 pack/case	16-2112121
2914 <b>4</b> 6	18" x 18"	5/pack, 40 pack/case	16-2118181
Rigid Pack			
McKesso <b>n #</b>	Description	Packaging	Mfr #
640985	18" x 18"	5/pack, 20 pack/case	16-2118182

## Trash Can with Plastic Liner entrust™ 32 Quart White Steel Step On Square

View Item



McKesson # 407341 UOM EA/1 Stock Stocked

Actual Product May Vary. See Product Details and Features for Specifications

#### **Product Features**

- entrust™ Waste Can
- Steel outer can and lid.
- Square
- 11-1/2D X 12-1/4W X 18-1/4H Inch
- 32-quart capacity (8 gallon).
- White enamel finish.
- Steel cans come with removable plastic inner liner and are available in enamel or stainless steel finishes.
- Step-on foot pedat opens lid.
- Perfect for exam rooms, restrooms, and front office areas.
- Meets OSHA standard on bloodborne pathogens.
- Hazardous waste label included.
- Not Made With Natural Rubber Latex.
- Packaged: Each

Manufacturer	McKesson Brand
Manufacturer #	81-35266
Brand	Entrust™
Application	Trash Can with Plastic Liner
Capacity	32 Quart
Color	White
Lid Opening Type	Step On
Lid Type	Self-Closing
Material	Steel
Size	11-1/2D X 12-1/4W X 18-1/4H Inch
Style	Square
UNSPSC Code	47121702
Latex Free Indicator	Not Made with Natural Rubber Latex

#### Questions? Contact Customer Service @ 866-625-2679

Login to SupplyManager Privacy Policy Terms of Use

More from McKesson Relay Health Primary Care Surgery Center

© 2016 McKesson Corporation

#### MCKESSON MEDICAL-SURGICAL, INC. - PROPRIETARY AND CONFIDENTIAL INFORMATION

Confidentiality Notice: This e-mail message, including any attachments, is for the sole use of the intended recipient(s) and may contain confidential and proprietary information. Any unauthorized review, use, disclosure or distribution is prohibited. If you are not the intended recipient, please contact the sender at <a href="SupplyManager@mckesson.com">SupplyManager@mckesson.com</a> and then permanently destroy all