QPatch 16X Screening Station by

Sophion Bioscience A/S Baltorpvej 154 DK-2750 Ballerup

Denmark

www.sophion.com

by sophion Vacuum: min. -700 mbar

Pressure: 4–8 bar 100–240V ~ 8,0–3.5A

50-60Hz

Fuses: T6.3AH/250Vac Made in Denmark

 \sim

This equipment complies with the requirement in Part 15 of FCC Rules. See manual. Protected by patents and patents pending. For patent information see http://www.sophion.dk/Technology/Patents